

Ultrasonic height and weight scale

MODEL : BYH01



SPECIFICATION

Function :

- >Height measuring function
- >Body fat measuring and analyzing function
Water/Muscle / Body bone / Calories display / BMI function

- > Material: ABS Plastic+ Aluminium alloy
- > Backlight LCD display: 102.5mm*67.5mm
- > Max weight: 5-200kg
- > Height measurement range: 80-195cm
- > Power : 4xAAA battery (included) or 6V DC

Packing details:

Product size: 44.8x34x 207.5cm

Color box size: 68.3 x37.5x22.5cm

Carton size: 68.3 x37.5x22.5cm

1 pcs/CTN

N.W.: 4.55 kg

G.W.: 5.93 kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel





Home



School

Ultrasonic Height Scale

MODEL : BYH02



(picture1)

(picture2)

SPECIFICATION

Function: Ultrasonic Height / Temperature function

- > Material: ABS Plastic
- > Backlight LCD display: 29*39mm
- > Max height: 500cm
- > Division: 0.5cm
- > Power: 4 x AAA battery (included)
- > Feature: Auto switch, auto zero, zero off, low battery / over-load indication

Packing details:

Product size: 50x59x 235mm

Color box size: 70x60x250mm

Carton size: 375x365x270cm

30 pcs/CTN

N.W.: 3.15kg

G.W.: 6.8kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel



Home



School

Height and Weight scale

MODEL : BYH03



SPECIFICATION

- Function:
- > Height measuring function
 - > Body fat measuring and analyzing function
Water/Muscle/Body bone/BMI/Calories display
 - > Voice function/Bluetooth function
 - > Connect with computer with RS232 port
-
- > Material: ABS glass+Aluminium alloy
 - > Blue backlit LCD display: 70*80mm
 - > Height measuring range: 0.80-1.95 m(Division:5mm)
 - > Weight measuring range: 0.3-200kg
 - > Power: 4*NO.7AAA batteres or 6V DC .With RS-232 port
 - > Feature: Step on to start,continuous weighing/Low battery / over-load indication/ Auto-zero

Packing details:

Product size: 35.5*49*210.4cm

Color box size: 123.5 x41x20.5cm

Carton size: 123.5 x41x20.5cm

1 PC/CTN

N.W.: 11.5kg

G.W.: 14.3kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel




Home



School

Height and weight scale with Printer

MODEL : BYH05



SPECIFICATION

Function:

- > Height measuring function/Body fat measuring and analyzing function
 - *Water/Muscle/Body bone/Calories display/BMI
- > Voice function
- > Connect with computer with RS232 port
- > Result printing function: With a thermal printer, press the print button to print the measurement result

- > Material: ABS plastic+ Aluminum
- > Blue backlit LCD display: 90*90mm
- > Height measuring range: 0.80-2.15 m(Division:2mm)
- > Weight measuring range: 5kg-200kg
- > Power: 6V DC power
- > Supply with RS-232 port
- > Feature: Step on to start /Low battery /over-load indication/ Auto-zero

Packing details:

Product size: 40*52.9*233cm

Color box size: 123.5 x55.8x22cm

Carton size: 123.5 x55.8x22cm

1 PC/CTN

N.W.: 18.5kg

G.W.: 22kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel




Home



School

MODEL : BYH07



SPECIFICATION

Function:

- > Height measuring function
- > Body fat measuring and analyzing function
Water/Muscle/Body bone/Calories display/BMI
- > Voice function
- > Result printing function: With a thermal printer, press the print button to print the measurement result

- > Material: glass+Aluminium alloy
- > Backlit LCD display size: 66*42mm
- > Height measuring range: 1.35-2.15 m (division of height :5mm)
- > Weight measuring range: 3kg-200kg
- > Power: 6V 3A DC power
- > Feature: Step on to start/Low battery / over-load indication
Auto-zero/Electricity saving function

Packing details:

Product size: 62*47*230cm

Color box size: 120*62*24.5cm

Carton size: 120*62*24.5cm

1 PC/CTN

N.W.: 20kg

G.W.: 22.60kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel



Home



School

Height and weight scale

MODEL : BYH09



SPECIFICATION

Function:

- > Height measuring function
- > Body fat measuring and analyzing function
Water/Muscle/Body bone/BMI
- > Calories display
- > Voice function

- > Material: glass+Aluminium alloy
- > Blue backlit LCD display size: 66*42.5mm
- > Height measuring range: 1.10-2.10 m(Division:5mm)
- > Weight measuring range: 3kg-200kg
- > Power: 6V 3A DC power
- > Feature: Step on to start,continuous weighing/Low battery / over-load indication
Auto-zero

Packing details:

Product size: 35.5*49*230.4cm

Color box size: 126.5*44.5*20cm

Carton size: 126.5*44.5*20cm

1 PC/CTN

N.W.: 11.6kg

G.W.: 13.20kg

Application:



Hospital



Beauty salon



GYM



Pharmacy



Supermarket



Hotel



Home



School

MODEL : BYHJ02



SPECIFICATION

Function: Height measuring function

> Material: ABS Plastic

> View window: 37*29mm

> Height measuring range: 0.70-2.10meter (division of height : 1mm)

Packing details:

Product size: 42*32*215cm

Carton size: 92*37*13cm

1 pc/CTN

N.W.: 3.00kg

G.W.: 2.60kg