



## COMPANY INTRODUCTION

Morntrip was established in 2015, which is the sales department of Xiantao Tongda Non-woven Products Co., Ltd. (since year of 1998). Our company located in Wuhan city, the biggest city in central China, Morntrip is a professional sales team for disposable non-woven products.

Morntrip's production bases are located in Xiantao city – the world famous and the biggest non-woven industry base in China. We have been in non-woven field for more than 20 years. At present, we have two production bases, one is for disposable clothing, one for disposable face mask and raw material. Two industry parks are more than 60,000 m<sup>2</sup>.

Morntrip produce non-woven raw material and non-woven finished products. Raw material including spun-bond PP, melt-blown PP and SMS. And finished non-woven products including disposable face mask, caps, gowns, coverall, shoe cover, bed sheet and etc. We can finish over 30\*40HQ containers of finished products and 400 tons of raw material every month. In 2017, our annual sales amount passed 10,000,000 USD.

As a manufacturer, we are certificated with national medical equipment production license, CE, ISO9001, ISO13485, NELSON and FDA. Our factory is one of the most important suppliers in the

year of 2002(SARS virus), year 2014(Ebola virus), and year 2015(MERS virus).

In the year of 2017, our new factory imported new automatic production lines from Japan for producing non-woven fabric, and we have built a few dust-free workshops for disposable face mask products. We are growing up years by years.

Now Morntrip's clients are from Europe, North America, Japan and Korea... and domestic market. We always have confidence in our products.

Morntrip always holding quality and honesty as the most important principles for our clients, which helps us growing up together with our customers. In the further, Morntrip will keep our spirit as always, pursuing innovation, providing best quality products and service for you!





# PRODUCTION WORKSHOPS 生产车间

原材料生产线 / Raw material production line



仓库 / Warehouses



生产车间 / Production Workshop





## QUALIFICATION HONOR 企业资质

We are qualified with CE, ISO 13485, FDA, NELSON and NIOSH certificates.



## CONTENTS

### 01. RAW MATERIAL

- Spun-bond Polypropylene
- Melt-blown Polypropylene
- SMS

### 02. MASK SERIES

#### Disposable Face Mask

- Earloop Face Mask
- Surgical Mask with Ties
- Active Carbon Face Mask
- Custom Printed Face Mask
- Face Mask with Shield

#### Disposable Dust Mask/Respirator

- Cone Shape Respirator (N95/FFP1/FFP2/FFP3)
- Respirator with Adjustable Straps
- Vertical Fold Respirator
- Flat Fold Respirator
- Nuisance Dust Mask
- Foldable Respirator with Valve

### 03. CAP SERIES

- Mob Cap
- Bouffant Cap
- Surgeon Cap with Ties
- Doctor Cap with Elastic
- Disposable Beard Cover

### 04. CLOTHING SERIES

#### Disposable Coverall

- PP/SMS Coverall
- Microporous Coverall
- Coverall With Taped Seams
- SMS Coverall With Bound Seams
- Microporous Coverall with SMS Back Panel

#### Disposable Gowns

- Isolation Gown
- Scrub Suit Set
- Surgical Gown
- Ultrasonic Heat Sealing Surgical Gown
- Reinforced Surgical Gown

#### Disposable Lab Coat

- Zipper Closure Lab Coat
- Snap Closure Lab Coat

#### Disposable Spa Beauty Items

- Kimono/Spa Robe/Bathrobe
- Disposable Bra
- Disposable Underwear
- Disposable Headband
- U Shape Pillow Cover

### 05. SHOE COVER SERIES

- PE/CPE/Non-woven Shoe Cover
- Anti-Slip Non-Woven Shoe Cover
- Anti-Static Non-Woven Shoe Cover
- Shoe Cover Coated with CPE Film

### 06. BED SET SERIES

- Disposable Bed Sheet/Cover
- Disposable Pillow Case
- Disposable Nursing Pad
- Disposable Pet Pad

### 07. OTHER PRODUCTS

- Restaurant Bib
- Disposable Oversleeves
- Disposable Gloves
- PU Foam/Silicone Earplugs



## Non-woven Raw Material



### Spun-bond Polypropylene

- Material: Polypropylene
  - Color: ☐ ☐ ☐ ☐ ☐
  - Fabric Width:  $\leq 330\text{cm}$
  - Packing: 100-500kg/roll
- Raw material for disposable face mask, caps, gowns, coverall, lab coat, shoe cover, bed sheet, non woven bags, etc.



### Melt-blown Polypropylene


- Material: Polypropylene
  - Color: ☐ ☐
  - Fabric Width:  $\leq 330\text{cm}$
  - Packing: 100-500kg/roll
- Raw material(as filter) for disposable face mask, dust mask, disposable respirators, etc.



### SMS

- Material: Polypropylene
  - Color: ☐ ☐ ☐
  - Fabric Width:  $\leq 330\text{cm}$
  - Packing: 100-500kg/roll
- Raw material for disposable face mask, caps, gowns, coverall, lab coat, shoe cover, bed sheet, non woven bags, etc.

## Disposable Face Mask



### 2 Ply Earloop Face Mask

- Model NO: MT-ME01
- Material: 2 ply spun-bond polypropylene
- Color: ☐ ☐ ☐ ☐ ☐
- Style: Earloop
- Elastic: 40D-70D round elastic
- Adult Size:  $17.5 \times 9.5\text{cm}$
- Teenage Size:  $14.5 \times 9.5\text{cm}$
- Kids Size:  $12 \times 7\text{cm}$
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.



### 3 Ply Earloop Face Mask

- Model NO: MT-ME02
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Color: ☐ ☐ ☐ ☐ ☐
- Style: Earloop
- Elastic: 40D-70D round elastic
- Adult Size:  $17.5 \times 9.5\text{cm}$
- Teenage Size:  $14.5 \times 9.5\text{cm}$
- Kids Size:  $12 \times 7\text{cm}$
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.




### 3 Ply Face Mask With Flat Earloop

- Model NO: MT-ME02F
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Color: ☐ ☐ ☐ ☐ ☐
- Style: Earloop
- Elastic: 2mm-4mm flat elastic
- Adult Size:  $17.5 \times 9.5\text{cm}$
- Teenage Size:  $14.5 \times 9.5\text{cm}$
- Kids Size:  $12 \times 7\text{cm}$
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.




## Disposable Face Mask



### 4 Ply Active Carbon Face Mask

- Model NO: MT-ME03
- Material: 4 ply non woven(SBPP+active carbon+MBPP+SBPP)
- Color: ■ ■
- Style: Earloop
- Elastic: Round/Flat Elastic
- Adult Size: 17.5 × 9.5cm
- Teenage Size: 14.5 × 9.5cm
- Kids Size: 12 × 7cm
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.



### 3 Ply Surgical Face Mask With Ties

- Model NO: MT-MT01
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Color: □ ■ ■
- Style: With ties
- Size: 17.5 × 9.5cm
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Mainly for medical use



### Custom Printed Face Mask

- Model NO: MT-MP01
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Printing Type: Full printing
- Style: Earloop
- Elastic: Round/Flat Elastic
- Adult Size: 17.5 × 9.5cm
- Teenage Size: 14.5 × 9.5cm
- Kids Size: 12 × 7cm
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.



### Custom Printed Face Mask

- Model NO: MT-MP02
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Printing Type: Registration printing
- Style: Earloop
- Elastic: Round/Flat Elastic
- Adult Size: 17.5 × 9.5cm
- Teenage Size: 14.5 × 9.5cm
- Kids Size: 12 × 7cm
- Standard Packing: 50pcs/box, 2000pcs/carton
- Application: Medical, daily use, household duties, cleaning, beauty, food industry and other industries.



### Procedure Face Masks with Shield

- Model NO: MT-ME02S
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Shield: Transparent plastic
- Color: □ ■ ■
- Style: Earloop
- Elastic: Round/Flat Elastic
- Size: Shield: 30 × 12.5cm, Mask: 17.5 × 9.5cm
- Standard Packing: 1pc/box, 50pcs/box, 1000pcs/carton
- Application: Mainly used in OR, during surgeries.



### Standard Surgical Masks with Shield

- Model NO: MT-MT01S
- Material: 3 ply non woven(SBPP+MBPP+SBPP)
- Shield: Transparent plastic
- Color: □ ■ ■
- Style: With ties
- Size: Shield: 30 × 12.5cm, Mask: 17.5 × 9.5cm
- Standard Packing: 1pc/box, 50pcs/box, 1000pcs/carton
- Application: Mainly used in OR, during surgeries.



### Plastic Face Mask

- Model NO: MT-MF01
- Material: Plastic
- Color: Transparent
- Style: Earloop
- Elastic: Round elastic
- Size: 14 × 6cm
- Standard Packing: 1pc/opp bag, 20pcs/box, 200pcs/carton
- Application: Mainly used in catering industry.



### 3D Fold Mask

- Model NO: MT-MZ01
- Material: Non-woven
- Color: □ ■ ■
- Style: Earloop
- Size: 12.5 × 11.5cm
- Standard Packing: 50pcs/box, 1000pcs/carton
- Application: Household duties, daily use, cleaning, beauty industry.



## Disposable Dust Mask/Respirator



### » Cone Shape Respirator

- Model NO: ZYB-11 / ZYB-11A
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: Needle-punched +filter+needle-punched nonwoven
- Color: ☐
- Head Straps: Nylon elastic
- Style: Ultrasonic welded head straps, aluminum nose clip
- Valve: With/without valve
- Standard Packing: 20pcs/box, 20boxes/carton
- Application: Medical, daily use, household duties, gardening, mining, sanding, drilling, cutting, plastering, pharmaceutical industry.



### » Cone Shape Respirator With Valve

- Model NO: ZYB-12 / ZYB-10(with active carbon)
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+needle-punched nonwoven
- Color: ☐ ☐ ☐
- Head Straps: Nylon elastic
- Style: Ultrasonic welded head straps, aluminum nose clip
- Valve: With/without valve
- Standard Packing: 20pcs/box, 20boxes/carton
- Application: Medical, daily use, household duties, gardening, mining, sanding, drilling, cutting, plastering, pharmaceutical industry.



### » Cone Shape Respirator

- Model NO: ZYB-17 / ZYB-17A
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+needle-punched nonwoven
- Color: ☐
- Head Straps: Nylon elastic
- Style: Ultrasonic welded head straps, aluminum nose clip
- Valve: With/without valve
- Standard Packing: 20pcs/box, 20boxes/carton
- Application: Medical, daily use, household duties, gardening, mining, sanding, drilling, cutting, plastering, pharmaceutical industry.



### » Cone Shape Respirator

- Model NO: ZYB-18A / ZYB-18
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+needle-punched nonwoven
- Color: ☐
- Head Straps: Rubber elastic
- Style: Stapled head straps, aluminum nose clip
- Valve: With/without valve
- Standard Packing: 20pcs/box, 20boxes/carton
- Application: Medical, daily use, household duties, gardening, mining, sanding, drilling, cutting, plastering, pharmaceutical industry.



### » FFP3 Cone Respirator With Adjustable Straps

- Model NO: ZYB-22 / ZYB-21
- Protective Level: FFP3
- Material: PP+filter+needle-punched nonwoven
- Color: ☐
- Head Straps: Wide nylon elastic
- Style: Adjustable head straps with four adjusters, full face seal contours, plastic nose clip
- Valve: With valve
- Standard Packing: 10pcs/box, 20boxes/carton
- Application: Medical, industry, mining, sanding, drilling, cutting, plastering, pharmaceutical industry, etc.




### » Vertical Fold Respirator

- Model NO: ZYE-01 / ZYE-02
- Protective Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+hot air cotton+SBPP, black foam nose cushion
- Color: ☐
- Head Straps: Nylon elastic
- Style: Vertical fold, welded head straps, plastic nose clip
- Valve: With/without valve
- Standard Packing: 25pcs/box, 20boxes/carton
- Application: Medical, daily use, industry, household duties, food industry, gardening, etc.



Disposable Dust Mask/Respirator



Flat Fold Respirator With Valve

- Model NO: ZYE-06 / ZYE-06C(with active carbon)
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+hot air cotton+SBPP, black foam nose cushion
- Color: ☐
- Head Straps: Nylon elastic
- Style: Flat fold, adjustable head straps, plastic nose clip
- Valve: With/without valve
- Standard Packing: 25pcs/box, 20boxes/carton
- Application: Medical, daily use, industry, household duties, food industry, gardening, etc.

Nuisance Dust Mask

- Model NO: ZYA-01
- Protection Level: Basic dust proof
- Material: Melt-blown polypropylene
- Color: ☐ ☐ ☐
- Head Straps: Nylon elastic
- Style: Stapled/welded head strap, aluminum nose clip
- Standard Packing: 50pcs/box, 20boxes/carton
- Application: Food industry, household duties, gardenin, etc.



Vertical Fold Respirator

- Model NO: ZYE-08 / ZYE-08C(with active carbon)
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+filter+hot air cotton+SBPP, black foam nose cushion
- Color: ☐
- Head Straps: Nylon elastic
- Style: Vertical fold, welded head straps, plastic nose clip
- Valve: With/without valve
- Standard Packing: 25pcs/box, 20boxes/carton
- Application: Medical, daily use, industry, household duties, food industry, gardening, etc.

Foldable Respirator With Valve

- Model NO: ZYE-11
- Protection Level: N95/FFP1/FFP2/FFP3
- Material: SBPP+hot air cotton+filter+hot air cotton+SBPP, black foam nose cushion
- Color: ☐
- Head Straps: Nylon elastic
- Style: Foldable, adjustable head straps
- Valve: With valve
- Standard Packing: 20pcs/box, 20boxes/carton
- Application: Medical, daily use, industry, household duties, food industry, gardening, etc.



## Disposable Caps



### Single Elastic Mob Cap

- Model NO: MT-CM01
- Material: PP
- Weight: 10-20gsm
- Color: ☐ ☒ ☒ ☒ ☒ ☒
- Size: 18"-24"
- Style: Single elastic
- Standard Packing: 100pcs/bag, 2000pcs/carton
- Application: Medical, household duties, cleaning, beauty & food industry, labs, catering industry.



### Double Elastic Mob Cap

- Model NO: MT-CM02
- Material: PP
- Weight: 10-20gsm
- Color: ☐ ☒ ☒ ☒ ☒ ☒
- Size: 18"-24"
- Style: Double elastic
- Standard Packing: 100pcs/bag, 2000pcs/carton
- Application: Medical, household duties, cleaning, beauty & food industry, labs, catering industry.



### Bouffant Cap

- Model NO: MT-CB01
- Material: PP/SMS
- Weight: 10-35gsm
- Color: ☐ ☒ ☒ ☒ ☒ ☒
- Size: 18"-24"
- Standard Packing: 100pcs/bag, 2000pcs/carton
- Application: Medical, household duties, cleaning, beauty & food industry, labs, catering industry.



### Surgeon Cap With Tie

- Model NO: MT-CS01
- Material: PP/SMS
- Weight: 20-40gsm
- Color: ☐ ☒ ☒ ☒
- Size: 13 x 64cm
- Standard Packing: 100pcs/bag, 2000pcs/carton
- Application: Mainly for surgeons



### Doctor Cap With Elastic

- Model NO: MT-CS02
- Material: PP/SMS
- Weight: 20-40gsm
- Color: ☐ ☒ ☒ ☒
- Size: 15 x 62cm
- Standard Packing: 100pcs/bag, 2000pcs/carton
- Application: Medical use, for doctors/nurses.



### Disposable Beard Cover

- Model NO: MT-CBC01
- Material: PP
- Weight: 10-20gsm
- Color: ☐ ☒ ☒ ☒
- Size: 10"-24"
- Seam: Heat welded
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Food industry, restaurant, catering industry.



## Disposable Coverall



### » Polypropylene Coverall

- Model NO: MT-CA02A
- Material: 20-65gsm PP
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Style: Hooded without boot, zipper closure without flap, elastic back waist, wrists and ankles.
- Size: S-5XL
- Stitch: Four threads overlock
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Cleanrooms, common dust proof, food industry, etc.



### » Polypropylene Coverall Coated With PE Film

- Model NO: MT-CA02B
- Material: 20-65gsm PP
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Style: Hooded without boot, zipper closure with flap, elastic wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Cleanrooms, dust proof, water proof, paint sprays, chemical protection, etc.



### » SMS Coverall With Shirt Collar

- Model NO: MT-CA01
- Material: 20-65gsm SMS
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Style: Shirt collar without boot, zipper closure with flap, elastic back waist, wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Cement manufacturing, grinding & polishing, light industrial cleaning, machine maintenance, etc.



### » SMS Coverall

- Model NO: MT-CA02C
- Material: 20-65gsm SMS
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Style: Hooded without boot, zipper closure without flap, elastic back waist, wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Cement manufacturing, grinding & polishing, light industrial cleaning, machine maintenance, etc.



## Disposable Coverall



### » Microporous Coverall

- Model NO: MT-CA02D
- Material: 40-65gsm microporous
- Color:
- Style: Hooded without boot, zipper closure with flap, elastic back waist, wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection and removal, site preparation, insulation, cement manufacturing, grinding & polishing, pharmaceutical processing & handling toxic powders, cleanrooms.



### » Microporous Coverall With Shoe Cover

- Model NO: MT-CA03
- Material: 50-70gsm microporous
- Color:
- Style: Hooded with boots, zipper closure with flap, elastic back waist, wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection and removal, site preparation, insulation, cement manufacturing, grinding & polishing, pharmaceutical processing & handling toxic powders, cleanrooms.





## Disposable Coverall



### ✧ Coverall With Taped Seams

- Model NO: MT-CA04
- Material: Microporous/SMS
- Weight: 50-65gsm
- Color: ☐ ☐ ☐ ☐ ☐
- Style: Hooded with boots, zipper closure, storm flap with adhesive tape, elastic wrists and back waist.
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 4&5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection, site preparation, insulation, cement manufacturing, grinding & polishing, industrial cleaning and maintenance, painting & coating, resin.



### ✧ SMS Coverall With Bound Seams

- Model NO: MT-CA05
- Material: SMS/Microporous
- Weight: 50-65gsm
- Color: ☐ ☐ ☐ ☐ ☐
- Style: 3-panel hood, zipper closure, storm flap with adhesive tape, elastic ankles, wrists and back waist.
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection and removal, site preparation, insulation, cement manufacturing, grinding & polishing, pharmaceutical processing & handling toxic powders, cleanrooms.





Disposable Coverall



Four Layer Aluminium Film Coverall

- Model NO: MT-CA07
- Material: PP+PE+aluminium+PE
- Weight: 80gsm
- Color: Silver
- Style: Hooded without boot, zipper closure with flap, elastic back waist, wrists and ankles
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 4&5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Emergency warm keeping, fire resistant, heat resistant.



Coverall With Reflective Tape

- Model NO: MT-CA08
- Material: Microporous/SMS
- Weight: 50-65gsm
- Color:
- Style: Hooded without boot, zipper closure with flap, elastic back waist, wrists and ankles
- Reflective Tape: Around chest, wrists and legs
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection and removal, site preparation, insulation, cement manufacturing, grinding & polishing, pharmaceutical processing & handling toxic powders, cleanrooms.



Microporous Coverall Combined With SMS Back Panel

- Model NO: MT-CA06
- Material: Microporous combined with SMS
- Weight: 50-65gsm
- Color:
- Style: Hooded without boot, zipper closure with flap, elastic back waist and ankles, knitted cuffs.
- Size: S-5XL
- Stitch: Four threads overlock
- Protection Level: Type 5&6
- Standard Packing: 1pc/bag, 50pcs/carton
- Application: Asbestos inspection and removal, site preparation, insulation, cement manufacturing, grinding & polishing, pharmaceutical processing & handling toxic powders, cleanrooms.





## Disposable Gowns



### » Elastic Cuff Isolation Gown

- Model NO: MT-IG01
- Material: PP/SMS/PP+PE/SMS+PE
- Weight: 20-40gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐
- Collar: Tie-on
- Cuffs: Elastic
- Waist: 2 ties
- Stitch: Three threads overlock
- Size: S-XXXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Low risk medical procedures, visiting patients, clean rooms, labs.



### » Impervious Knitted Cuffs Isolation Gown

- Model NO: MT-IG01VK
- Material: PP/SMS/PP+PE/SMS+PE
- Weight: 20-40gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐
- Collar: Velcro
- Cuffs: Knitted
- Waist: 2 ties
- Stitch: Three threads overlock
- Size: S-XXXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Low risk medical procedures, visiting patients, clean rooms, labs.



### » Short Sleeves Isolation Gown

- Model NO: MT-IG02
- Material: PP/SMS/PP+PE/SMS+PE
- Weight: 20-40gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐
- Collar: Tie-on
- Waist: 2 ties
- Stitch: Three thread overlock
- Size: M-XXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Mainly for patients, health examination.



### » No Sleeve Isolation Gown

- Model NO: MT-IG03
- Material: PP/SMS/PP+PE/SMS+PE
- Weight: 20-40gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐
- Collar: Tie-on
- Waist: 2 ties
- Stitch: Three thread overlock
- Size: S-XXXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Mainly for patients, health examination.



## Disposable Gowns



### » Scrub Suit Set

- Material: PP/SMS
- Weight: 25-50gsm
- Color:
- Collar: V or round collar
- Sleeves: Short sleeves
- Waist: With tie or elastic
- Stitch: Three thread overlock
- Size: S-XXXL
- Standard Package: 10sets/bag, 50sets/carton
- Application: Mainly used by doctors and nurses.



### » Surgical Gown

- Model NO: MT-SG01
- Material: SMS/SMMS
- Weight: 30-50gsm
- Color:
- Collar: Tie-on or velcro
- Cuffs: 7-10cm knitted
- Waist: 4 ties
- Stitch: Three thread overlock
- Size: S-XXXL
- Standard Packing: EO gas sterilized with paper-plastic bag packing.
- Application: Mainly wore by surgeons/nurses during surgeries.



### » Ultrasonic Heat Sealing Surgical Gown

- Model NO: MT-SG02
- Material: SMS/SMMS
- Weight: 35-50gsm
- Color:
- Collar: Velcro
- Cuffs: 7-10cm knitted
- Waist: 4 ties
- Seam: Ultrasonic heat sealing
- Size: S-XXXL
- Standard Packing: EO gas sterilized with paper-plastic bag packing.
- Application: Mainly wore by surgeons/nurses during surgeries.




## Disposable Gowns



### Reinforced Surgical Gown

- Model NO: MT-SG03
- Material: SMS/SMMS
- Weight: 35-50gsm
- Color: ■ ■
- Reinforced Part: Chest and sleeves
- Collar: Velcro
- Cuffs: 7-10cm knitted
- Waist: 4 ties
- Seam: Ultrasonic heat sealing
- Size: S-XXXL
- Standard Packing: Eo gas sterilized with paper-plastic bag packing.
- Application: Mainly wore by surgeons/nurses during surgeries.



### Standard Thumb CPE Gown

- Material: CPE
- Weight: 30-50g/pc
- Size: S-XXL
- Style: Long sleeve, opening cuff, with thumb loop, with ties at waist
- Color: □ ■
- Cuff: Opening cuff, with thumb loop
- Waist: 2 ties
- Stitch: Three thread overlock
- Standard Packing: 1pc/bag, 10pcs/middle bag, 100pcs/carton
- Application: Food processing, cleaning, beauty salon, medical & chemical tests, painting, medical use, butchery, etc.





## Disposable Lab Coat



### Zipper Closure Lab Coat

- Model NO: MT-LCZ01K
- Material: PP/SMS/PP+PE/SMS+PE
- Color: ☐ ☒ ☐ ☐ ☐ ☐
- Weight: 20-50gsm
- Closure: Zipper
- Collar: Knitted
- Cuffs: Elastic
- Pocket: Three pocket
- Stitch: Three thread overlock
- Size: M/L/XL/XXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Medical institutions, labs, clean room.



### Snap Closure Lab Coat

- Model NO: MT-LCP02K
- Material: PP/SMS/PP+PE/SMS+PE
- Color: ☐ ☒ ☐ ☐ ☐ ☐
- Weight: 20-50gsm
- Closure: Snap
- Collar: Knitted
- Cuffs: Knitted
- Pocket: None
- Stitch: Three thread overlock
- Size: M/L/XL/XXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Medical institutions, labs, clean room.



### Velcro Closure Lab Coat

- Model NO: MT-LCV01
- Material: PP/SMS/PP+PE/SMS+PE
- Color: ☐ ☒ ☐ ☐ ☐ ☐
- Weight: 20-50gsm
- Closure: Velcro
- Collar: Shirt collar
- Cuffs: Elastic
- Pocket: None
- Stitch: Three thread overlock
- Size: M/L/XL/XXL
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Medical institutions, labs, clean room.



## Disposable Spa Beauty Items



### Long Sleeve Kimono/Spa Robe

- Model NO: MT-KL01
- Material: PP/SMS/Spunlace
- Weight: 20-50gsm
- Style: Long sleeves
- Waist: Ties on
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Stitch: Three thread overlock
- Size: Customized
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Hotel, beauty salon, SPA, sauna.



### Short Sleeve Kimono/Spa Robe

- Model NO: MT-KS01
- Material: PP/SMS/Spunlace
- Weight: 20-50gsm
- Style: Short sleeves
- Waist: Ties on
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Stitch: Three thread overlock
- Size: Customized
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Hotel, beauty salon, SPA, sauna.



### Disposable Bathrobe

- Model NO: MT-KS02
- Material: PP/SMS
- Weight: 35-50gsm
- Style: With elastic
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Size: Customized
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: Hotel, beauty salon, SPA, sauna, etc.



### Disposable Bra

- Material: PP/SMS
- Weight: 30-45gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Size: Customized size
- Stitch: Three thread overlock
- Style: One piece with elastic shoulder ties
- Standard Packing: 1pc/PE bag, 5pcs/opp bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon, hospital.



### Disposable Underwear

- Material: PP/SMS
- Weight: 30-45gsm
- Color: ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐
- Size: Customized size
- Stitch: Three thread overlock
- Style: Triangle
- Standard Packing: 1pc/PE bag, 5pc/opp bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon, hospital, traveling, etc.



Disposable Spa Beauty Items



» Disposable Underwear with Printing Design

- Material: PP/Spunlace
- Weight: 25-60gsm
- Color: printed
- Size: Customized size
- Stitch: Four thread overlock
- Style: Triangle
- Standard Packing: 1pc/PE bag, 5pc/opp bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon, hospital, traveling, etc.



» Cotton Disposable Underwear

- Material: Cotton
- Weight: 90gsm
- Color:
- Size: Customized size
- Stitch: Five thread overlock
- Style: Triangle
- Standard Packing: 5pcs/PE bag, 600pcs/carton
- Application: SPA, sauna, beauty salon, hospital, traveling, etc.



» Disposable G-string

- Material: PP/SMS/Spunlace/Cotton
- Weight: 30-45gsm
- Color:
- Size: Customized size
- Stitch: Four thread overlock
- Style: G-string
- Standard Packing: 1piece/PE bag, 5pc/opp bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon, hospital.



» Disposable Shorts Pants/Boxer

- Material: PP/SMS
- Weight: 30-45gsm
- Color:
- Size: Customized size
- Stitch: Four thread overlock
- Style: Boxer shorts
- Standard Packing: 10pcs/bag, 200pcs/carton
- Application: SPA, sauna, beauty salon, hospital.



» Disposable Headband

- Model NO: MT-HB01
- Material: 2 ply PP
- Weight: 10gsm
- Style: 3-4 elastic
- Color:
- Size: 35 x 11cm
- Standard Packing: Individual packing, 100pcs/bag, 1000pcs/carton
- Application: Hotel, beauty salon, SPA, sauna, etc.



» Disposable Towel

- Material: Wood-pulp spunlace/paper
- Weight: 40-70gsm
- Color:
- Size: 40 x 60cm or customized size
- Standard Packing: 100pcs/bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon, hospital.



» U Shape Pillow Cover

- Material: PP
- Weight: 25-35gsm
- Color:
- Size: 36 x 33cm
- Standard Packing: 50pcs/bag, 500pcs/carton
- Application: SPA, sauna, beauty salon.



» Y Shape Face Rest Cover

- Material: Spunlace
- Weight: 45gsm
- Color:
- Size: 43 x 34cm
- Standard Packing: 100pcs/bag, 1000pcs/carton
- Application: SPA, sauna, beauty salon.



## Shoe Cover Series



### PE Shoe Cover

- Model NO: MT-SCM01
- Material: PE(polyethylene)
- Weight: 2-5g/pc
- Color: ■ Transparent
- Size: 38 × 16cm - 42 × 17cm
- Type: Machine made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.

### CPE Shoe Cover

- Model NO: MT-SCM01C
- Material: CPE(Polyvinyl chloride)
- Weight: 2.5-5g/pc
- Color: ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Machine made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.

### Machine Made Non-woven Shoe Cover

- Model NO: MT-SCM02
- Material: PP/SMS
- Weight: 20-40gsm
- Color: ■ ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Machine made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.

### Machine Made Non-woven Anti-slip Shoe Cover

- Model NO: MT-SCM03
- Material: PP/SMS with anti-slip grain
- Weight: 20-40gsm
- Color: ■ ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Machine made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.



### Hand Made Non-woven Shoe Cover

- Model NO: MT-SCH01
- Material: PP/SMS
- Weight: 20-40gsm
- Color: ■ ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Hand made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.

### Hand Made Non-woven Anti-slip Shoe Cover

- Model NO: MT-SCH02
- Material: PP/SMS with anti-slip grain
- Weight: 20-40gsm
- Color: ■ ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Hand made
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Hospital, clean rooms, labs, household, industry.

### Anti-static Non-woven Shoe Cover

- Model NO: MT-SCH03
- Material: PP/SMS with ESD tape
- Weight: 20-40gsm
- Color: ■ ■
- Size: 38 × 16cm - 42 × 17cm
- Type: Antistatics
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Labs, electricity industry or other industries.

### PP/SMS Shoe Cover Coated with CPE Film

- Model NO: MT-SCM02C
- Material: PP/SMS + CPE
- Weight: 20-40gsm
- Color: ■ ■ ■
- Size: 38 × 16cm - 42 × 18cm
- Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
- Application: Labs, electricity industry or other industries.



Disposable Bed Sets



Disposable Bed Sheet

- Material: PP/SMS/PP+PE/SMS+PE/spunlace
- Weight: 20-45gsm
- Color: ☐ ☒ ☒
- Size: Customized size
- Style: Flat
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: SPA, beauty, hospital.



Disposable Bed Cover

- Material: PP/SMS/PP+PE/SMS+PE/spunlace
- Weight: 20-45gsm
- Color: ☐ ☒ ☒
- Size: Customized size
- Style: Elastic on 2 ends, on 4 corners or all around
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: SPA, beauty, hospital.




Disposable Pillow Case

- Material: PP/SMS/PP+PE/SMS+PE/spunlace
- Weight: 20-45gsm
- Color: ☐ ☒ ☒
- Size: Customized size
- Closure: Open/flap/zipper/velcro
- Standard Packing: 10pcs/bag, 100pcs/carton
- Application: SPA, beauty, hospital.


Disposable Nursing Pad

- Material: Non-woven top sheet+fluff pulp & sanitary paper & SAP+PE film
- Size: 60 x 60cm, 60 x 90cm
- Color: ☐ ☒ ☒
- Weight: 60-110g/pc
- Standard Packing: 10-15pcs/bag, 12-20bags/carton
- Application: Nursing home, hospital, home supplies.



Disposable Pet Pad

- Material: Non-woven top sheet+fluff pulp & sanitary paper & SAP +PE film
- Size: 60 x 45cm, 60 x 60cm
- Weight: 40-90g/pc
- Color: ☐ ☒ ☒
- Standard Packing: 10-20pcs/bag, 10-20bags/carton
- Application: Nursing home, hospital, home supplies.





Other Products ➡



- » Plastic Restaurant Bib
- Model NO: MT-BRT01
  - Material: Plastic
  - Thickness: 0.04mm
  - Color: ☐
  - Logo: Customized
  - Size: 40 × 56cm with a 6cm pocket
  - Style: Tie on
  - Standard Packing: 100pcs/bag, 2500pcs/carton
  - Application: Catering industry, restaurant, hotel



- » Paper Plastic Bib
- Model NO: MT-BRT03
  - Material: Paper-plastic
  - Weight: 50-80gsm
  - Color: ☐
  - Logo: Customized
  - Size: Children: 39 × 57cm, Adults: 43 × 70cm
  - Style: Tie on
  - Standard Packing: 100pcs/bag, 1000pcs/carton
  - Application: Catering industry, restaurant, hotel.



- » PE/CPE Oversleeves
- Model NO: MT-SCP01
  - Material: PE/CPE
  - Size: 20 × 40cm, 20 × 44cm
  - Color: ☐ ☐ ☐ ☐ Transparent
  - Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
  - Application: Hospital, clean room, food industry, labs, restaurants, cleaning.



- » Non-woven Oversleeves
- Model NO: MT-SCN01
  - Material: PP/SMS
  - Size: 20 × 40cm, 20 × 44cm
  - Color: ☐ ☐ ☐ ☐
  - Standard Packing: 10pcs/roll, 10rolls/bag, 2000pcs/carton
  - Application: Hospital, clean room, food industry, labs, restaurants, cleaning.


- » SPA Liner
- Material: LDPE
  - Thickness: 0.012mm or 0.013mm
  - Size: S=105cm, M=112cm, L=118cm, XL=124cm
  - Color: ☐ ☐ ☐ Transparent
  - Standard Packing: 500pcs/carton
  - Application: SPA, beauty salon, pedicure.



- » Disposable Gloves
- Material: HDPE/PVC/Nitrile/Latex
  - Weight: Customized
  - Size: S/M/L/XL/XXL
  - Standard Packing: 100pcs/box, 1000pcs/carton
  - Application: Medical use, daily use, food industry, restaurants.



- » PU Foam Earplugs
- Material: PU foam
  - Color: ☐ ☐ ☐ ☐
  - Shape: Bullet/cylinder
  - Style: With/Without cord
  - Cord Material: PVC/Polyester
  - Cord Length: 70cm
  - Rebounded Time: 45-60seconds
  - Noise Reducing: 30-35dB
  - Size: 12 × 24mm, 13 × 24mm, 14 × 24mm, 13 × 18mm
  - Standard Packing: 1pair/bag, 100-200bags/box, 10-20boxes/carton
  - Application: Noise protection



- » Silicone Earplugs
- Material: Silicone
  - Color: ☐ ☐ ☐ ☐
  - Shape: Bullet
  - Style: With/Without cord
  - Cord Material: PVC/Polyester
  - Cord Length: 70cm
  - Noise Reducing: 29-33dB
  - Size: 11.5 × 27.5mm, 12.5 × 25.5mm, 12.5 × 31.5mm, 14.5 × 28.5mm, 15 × 31.5mm
  - Standard Packing: 1pair/bag, 100-200bags/box, 10-20boxes/carton
  - Application: Noise protection

