






In the field of laboratory products, Joanlab is a better choice for distributors and users. Joanlab optimizes the user experience by incorporating customer desired innovation and safety.

High quality product is the hallmark of Joanlab. We are able to provide low cost products to our customers without compromising our high quality and safety. In the past nine years we have incorporated innovative automation in our manufacturing processes and pass on the savings to our customers. We design, modify and manufacture ergonomically, user-friendly and environmentally-friendly products. Our easy to use and reliable products allow our customers to focus on their research and be more productive.

Joanlab started as a small company and overcome all the challenges to become a strong established laboratory products company. Joanlab has an experienced and collaborative R&D, manufacturing and marketing & sales teams that are committed to providing high quality and safe products for our customers. Thanks to our valuable employees, partners, and customers; Joanlab develops at more rapid speed and offer better services for more customers.





> Laboratory Glassware

1-36



> Laboratory Consumables

37-58



> Laboratory Equipment

59-90



> Laboratory Metal Ware

91-107

Laboratory Metal Ware





THREE FINGER CLAMP ALUMINUM SHAFT (BIG)

Item.	Name	Specification	Size
CJ170-10	Three finger clamp aluminum shaft (Big)	Chrome finish red*blue	Clamp70mm



THREE FINGER CLAMP COPPER SHAFT (BIG)

Item.	Name	Specification	Size
CJ170-18	Three finger clamp copper shaft (Big)	Chrome finish red*blue	Clamp70mm



THREE FINGER CLAMP(BIG)

Item.	Name	Specification	Size
CJ170-11	Three finger clamp(big)	Paint blue	Clamp70mm



THREE FINGER CLAMP (SMALL)

Item.	Name	Specification	Size
CJ170-20	Three finger clamp (small)	Chrome finish	Clamp25mm



THREE FINGER TURNING CLAMP(BIG)

Item.	Name	Specification	Size
SW300-40	Three finger turning clamp(Big)	Chrome finish with middle screw	Clamp70mm
SW300-45	Three finger turning clamp(Big)	Paint with middle screw, blue	Clamp70mm



THREE FINGER TURNING CLAMP(SMALL)

Item.	Name	Specification	Size
SW300-20	Three finger turning clamp(small)	Chrome finish with middle screw	Clamp25mm



THREE FINGER DUAL ADJUSTABLE CLAMP (LARGE)

Item.	Name	Specification	Size
FW200-70	Three finger dual adjustable clamp (large)	Double adjusting straight handle chrome finish	Clamp70mm
FW200-75	Three finger dual adjustable clamp (large)	Double adjusting straight handle, paint blue	Clamp70mm



THREE FINGER DUAL ADJUSTABLE CLAMP (SMALL)

Item.	Name	Specification	Size
FW200-30	Three finger dual adjustable clamp (small)	Double adjusting straight handle chrome finish	Clamp55mm
FW200-35	Three finger dual adjustable clamp (small)	Double adjusting straight handle, paint blue	Clamp55mm



THREE FINGER DUAL ADJUSTMENT TO MULTI CLAMP (LARGE)

Item.	Name	Specification	Size
SW200-70	Three finger dual adjustment to multi clamp (large)	Double adjustment chrome	Clamp70mm
SW200-75	Three finger dual adjustment to multi clamp (large)	Double adjusting straight handle, paint blue	Clamp70mm



THREE FINGER DUAL ADJUSTMENT TO MULTI USE CLIP (SMALL)

Item.	Name	Specification	Size
SW200-70	Three finger dual adjustment to multi use clip (small)	Double adjustment chrome	Clamp55mm
SW200-75	Three finger dual adjustment to multi use clip (small)	Double adjusting straight handle, paint blue	Clamp55mm



ENGLISH FOUR FINGER CLAMP

Item.	Name	Specification	Size
CJ190-10	English four finger clamp	chrome finish with rod	Clamp80mm
CJ190-11	English four finger clamp	blue paint with rod	Clamp80mm



EUROPEAN FOUR FINGER CLAMP

Item.	Name	Specification	Size
CJ200-10	European four finger clamp	chrome finish with rod	Clamp80mm
CJ200-11	European four finger clamp	blue paint with rod	Clamp80mm



ENGLISH FOUR FINGER TURNING CLAMP

Item.	Name	Specification	Size
SW300-80	English four finger turning clamp	chrome finish with middle screw	Clamp80mm
SW300-85	English four finger turning clamp	paint with middle screw, blue	Clamp80mm



EUROPEAN FOUR FINGER TURNING CLAMP

Item.	Name	Specification	Size
SW300-90	European four finger turning clamp	chrome finish with middle screw	Clamp90mm
SW300-95	European four finger turning clamp	paint with middle screw, blue	Clamp90mm



UNIVERSAL CLAMP 4-FINGERS

Item.	Name	Specification	Size
SJ200-10	Four finger turning clamp with rod	chrome finish	Clamp80mm
SJ200-11	Four finger turning clamp with rod	Paint	Clamp80mm



CHAIN CLAMP

Item.	Name	Specification	Size
CC400-145	Chain clip	Clamping diameter 40-145mm	clamp40mm-145mm



TWO FINGER FLASK CLAMP

Item.	Name	Specification	Size
FW100-45	Two finger flask clamp (Middle)	Middle chrome finish	Clamp12mm-45mm
FW100-25	Two finger flask clamp (Small)	Small chrome finish	Clamp7mm-25mm



TWO FINGER TURNING FLASK CLAMP

Item.	Name	Specification	Size
SW100-45	Two finger turning flask clamp (Middle)	Middle chrome finish	Clamp12mm-45mm
SW100-25	Two finger turning flask clamp (Small)	Small chrome finish	Clamp7mm-25mm



TWO FINGER DUAL ADJUSTABLE CLAMP

Item.	Name	Specification	Size
FT260-65	Two finger dual adjustable clamp (large)	Chrome finish	Clamp65mm
FT260-70	Two finger dual adjustable clamp (large)	Paint blue	Clamp65mm
FT260-40	Two finger dual adjustable clamp (small)	Chrome finish	Clamp40mm
FT260-45	Two finger dual adjustable clamp (small)	Paint blue	Clamp40mm



TWO FINGER DUAL ADJUSTMENT MULTI CLAMP

Item.	Name	Specification	Size
ST260-65	Two finger dual adjustment multi clamp (large)	Chrome finish	Clamp65mm
ST260-70	Two finger dual adjustment multi clamp (large)	Paint blue	Clamp65mm
ST260-40	Two finger dual adjustment multi clamp (small)	Chrome finish	Clamp40mm
ST260-45	Two finger dual adjustment multi clamp (small)	Paint blue	Clamp40mm



TWO FINGER DOUBLE HEAD CLAMP

Item.	Name	Specification	Size
TT600-10	Two finger double head clamp	Medium+Small	one side45mm, one side25mm

FLAST CLAMP

Item.	Name	Specification	Size
CJ204-28	Flast clamp	Alloy white paint	Clamp45mm
CJ204-30	Flast clamp	Spray-paint	Clamp55mm
CJ204-40	Flast clamp	Chrome finish	Clamp55mm



CJ204-28



CJ204-30



CJ204-40

UNIVERSAL CLAMP


Item.	Name	Specification	Size
CJ170-200	Iron universal clamp	Medium	Medium
CJ170-220	Iron universal clamp	Big	Large
CJ170-280	Iron universal clamp	XL	Outsize
CJ180-200	Copper universal clamp	Medium	Medium
CJ180-220	Copper universal clamp	Big	Large
CJ180-280	Copper universal clamp	XL	Outsize
CJ190-200	Stainless steel universal clamp	Medium	Medium
CJ190-220	Stainless steel universal clamp	Big	Large
CJ190-280	Stainless steel universal clamp	XL	Outsize




Iron universal clamp



Copper universal clamp



Stainless steel universal clamp



CULTURE DISH, EVAPORATING DISH TONGS

Item.	Name	Specification	Size
TB400-10	Culture dish, evaporating dish pliers	240mm	Length240mm
TB400-20	Culture dish, evaporating dish pliers	300mm	Length300mm




BEAKER TONGS

Item.	Name	Specification	Size
TB500-10	Beaker pliers	250mm	Length250mm
TB500-20	Beaker pliers	280mm	Length280mm



BEAKER TONGS (CHROME PLATED)

Item.	Name	Specification	Size
TB600-10	Beaker clamp (chrome plated) black rubber sleeve	260mm	Length260mm
CJ205-30	Beaker clamp (chrome plated) dip	260mm	Length260mm



CRUCIBLE TONGS

Item.	Name	Specification	Size
TA500-200	Crucible tongs	Chrome finish 200mm 8"	Length200mm
TA500-250	Crucible tongs	Chrome finish 250mm 10"	Length250mm
TA500-300	Crucible tongs	Chrome finish 300mm 12"	Length300mm
TA500-350	Crucible tongs	Chrome finish 350mm 14"	Length350mm
TA500-400	Crucible tongs	Chrome finish 400mm 16"	Length400mm
TA500-450	Crucible tongs	Chrome finish 450mm 18"	Length450mm
TA500-500	Crucible tongs	Chrome finish 500mm 20"	Length500mm
TA500-600	Crucible tongs	Chrome finish 600mm	Length600mm
TA600-150	Crucible tongs	Stainless steel 150mm	Length150mm
TA600-200	Crucible tongs	Stainless steel 200mm	Length200mm
TA600-250	Crucible tongs	Stainless steel 250mm	Length250mm
TA600-300	Crucible tongs	Stainless steel 300mm	Length300mm
TA600-350	Crucible tongs	Stainless steel 350mm	Length350mm
TA600-400	Crucible tongs	Stainless steel 400mm	Length400mm
TA600-450	Crucible tongs	Stainless steel 450mm	Length450mm
TA600-500	Crucible tongs	Stainless steel 500mm	Length500mm
TA600-600	Crucible tongs	Stainless steel 600mm	Length600mm



TEST TUBE CLAMP

Item.	Name	Specification	Size
TH800-10	Test tube clamp	Without handle	Clamp30mm
TH800-20	Test tube clamp	With handle	Clamp30mm
TH810-10	Wooden test tube clamp	Standard	Clamp20mm



TEST TUBE CLAMP

Item.	Name	Specification	Size
FD100-25	Alloy test tube clamp	Alloy chrome	Clamp25mm



THERMOMETER CLAMP (STRAIGHT ROD)

Item.	Name	Specification	Size
SD100-25	Alloy turnig test tube clamp	Alloy chrome	Clamp25mm



THERMOMETER CLAMP (STRAIGHT ROD)

Item.	Name	Specification	Size
ST200-15	Thermometer clamp (straight rod)	304 stainless steel	Extension rod 10mm



THERMOMETER CLAMP (STRAIGHT ROD)

Item.	Name	Specification	Size
ST120-15	Turning temperature thermometer clamp	304 stainless steel	Extension rod 10mm



GERMAN BOSS HEAD CLIP XL

Item.	Name	Specification	Size
CB200-08	German boss head clip XL	White paint	Clamp19mm



GERMAN BOSS HEAD CLIP

Item.	Name	Specification	Size
CB200-10	German boss head clip	Black paint	Clamp16mm



BOSS HEAD CLIP

Item.	Name	Specification	Size
CB202-13	Boss head clip	White, 13mm	Clamp13mm
CB202-16	Boss head clip	White, 16mm	Clamp16mm
CB202-21	Boss head clip	White, 21mm	Clamp21mm
CB202-25	Boss head clip XL	White, 25mm	Clamp25mm



BOSS HEAD CLIP

Item.	Name	Specification	Size
CB203-13	Boss head clip	Black, 13mm	Clamp13mm
CB203-16	Boss head clip	Black, 16mm	Clamp16mm
CB203-21	Boss head clip	Black, 21mm	Clamp21mm
CB203-25	Boss head clip XL	Balck , 25mm	Clamp25mm

**BOSS HEAD CLIP**

Item.	Name	Specification	Size
CB201-13	Boss head clip	Chrome finish, 13mm	Clamp13mm
CB201-16	Boss head clip	Chrome finish, 16mm	Clamp16mm
CB201-21	Boss head clip	Chrome finish, 21mm	Clamp21mm
CB201-25	Boss head clip XL	Chrome finish, 25mm	Clamp25mm

**BRITISH BOSS HEAD CLIP**

Item.	Name	Specification	Size
CB200-12	British boss head clip	Paint	Clamp15mm
CB200-15	British boss head clip	Chrome finish	Clamp15mm

**COPPER BOSS HEAD CLIP**

Item.	Name	Specification	Size
CB200-16	Copper boss head clip	Medium size	--

**TURNING BOSS HEAD**

Item.	Name	Specification	Size
BB180-16	Turning boss head	16mm	Clamp16mm

**SQUARE BOSS HEAD**

Item.	Name	Specification	Size
CB280-10	Square boss head	10*10	Clamp10*10mm
CB280-12	Square boss head	12*12	Clamp12*12mm

**BALL SHAPE BOSS HEAD**

Item.	Name	Specification	Size
CB270-13	Ball shape boss head	--	Clamp7-13mm

**MULTI USE THREE DIRECTION BOSS HEAD CLIP**

Item.	Name	Specification	Size
CB290-10	Multi use three direction boss head clip	Paint	Clamp12mm
CB290-20	Multi use three direction boss head clip	Chrome finish	Clamp12mm



CB180-10/CB180-12



CB170-10



CB170-20



CB170-30/CB170-31

BURETTE CLAMP

Item.	Name	Specification
CB180-10	German titration clip	double side plastic
CB180-12	German titration clip	single side plastic
CB170-10	German titration clip	double side alu
CB170-20	German titration clip	single side alu
CB170-30	German titration clip	single side chrome finish
CB170-31	German titration clip	single side chrome finish

**CORK BORERS WITH HANDLE (IRON PLATING)**


Item.	Name	Specification	Size
CB570-4	Punch (iron plating)	4*1	--
CB570-6	Punch (iron plating)	6*1	--
CB570-8	Punch (iron plating)	8*1	--

**CORK BORERS WITH HANDLE (COOPER TEETH)**

Item.	Name	Specification	Size
CB580-03	Punch (copper teeth)	5-7.5mm	5-7.5mm
CB580-06	Punch (copper teeth)	5-12.5mm	5-12.5mm
CB580-12	Punch (copper teeth)	5-18.75mm	5-18.75mm
CB580-18	Punch (copper teeth)	5-26.25mm	5-26.25mm

**IRON STAND MIDIEM SIZE**

Item.	Name	Size(mm)	Rod(mm)
AA330-140	Iron stang midiem size	200*140	600*10

**IRON STAND BIG SIZE**

Item.	Name	Size(mm)	Rod(mm)
AA330-180	Iron stang big size	300*180	700*10

**304 STAINLESS STAND**

Item.	Name	Size(mm)	Rod(mm)
AA330-140A	304 stainless stand	200*140	580*10



BRITISH IRON STAND

Item.	Name	Size(mm)	Rod(mm)
AA390-260	British iron stand	200*125	700*10
AA390-270	British iron stand	250*160	700*10
AA390-280	British iron stand	315*200	700*10
AA390-310	British iron stand	400*250	700*10



BRITISH PUNCH STAND (IRON PAINTING)

Item.	Name	Size(mm)	Rod(mm)
AA380-175	British punch stand (iron painting)	175*100 paint	700*10
AA380-225	British punch stand (iron painting)	225*125 paint	700*10
AA380-230	British punch stand (iron painting)	230*130 paint	700*10
AA380-250	British punch stand (iron painting)	250*150 paint	700*10
AA380-300	British punch stand (iron painting)	300*150 paint	




BRITISH PUNCH STAND(STAINLESS STEEL)


Item.	Name	Size(mm)	Rod(mm)
AA400-175	British punch stand (stainless steel)	175*100 stainless steel	700*10
AA400-225	British punch stand (stainless steel)	225*125 stainless steel	700*10
AA400-230	British punch stand (stainless steel)	230*130 stainless steel	700*10
AA400-250	British punch stand (stainless steel)	250*150 stainless steel	700*10
AA400-300	British punch stand (stainless steel)	300*150 stainless steel	700*10

STAND

Item.	Name	Size(mm)	Rod(mm)
AA300-150	Marble stand	300*150	650*10
AA360-190	Glass stand	300*150Toughened glass, iron frame	650*10
CA300-150	Ceramic stand	300*150	650*10



AA300-150



AA360-190



CA300-150



THREE BAR TYPE TABLE

Item.	Name	Size(mm)	Rod(mm)
AA100-010	Three bar type table	500*250	800*10
AA200-010	Three bar type table	400*200	700*10

SQUARE BRACKET (COMMON TYPE)

Item.	Name	Specification	Iron Coated Size(mm)	Rod(mm)
SP200-140	Square bracket (common type)	200*140mm	60mm/80mm/100mm	600*10
SP200-140S	Square bracket (common type)	200*140mm	60mm/80mm/100mm	580*10



SP200-140



SP200-140S



THREE LEGGED TABLE

Item.	Name	Size(mm)	Rod(mm)
AA320-85	Three legged tabl	85	--
AA320-110	Three legged tabl	110	--
AA320-135	Three legged tabl	135	--
AA320-175	Three legged tabl	175	--



A TYPE TABLE

Item.	Name	Specification	Rod(mm)
AA420-100	A type table	Small	610*10
AA420-150	A type table	Medium	610*10
AA420-200	A type table	Big	610*10



ALUMINIUM OXIDE LAB-JACKS

Item.	Name	Specification	Size(mm)
LJ450-100	Aluminium oxide lab-jacks	--	100*100
LJ450-150	Aluminium oxide lab-jacks	--	150*150
LJ450-200	Aluminium oxide lab-jacks	--	200*200
LJ450-250	Aluminium oxide lab-jacks	--	250*250



LJ480-150



LJ480-200



LJ550-100



LJ550-150



LJ550-200


STAINLESS STEEL LAB-JACKS

Item.	Name	Specification	Size(mm)
LJ480-150	430 stainless steel lab-jacks	430 stainless steel	150*150
LJ480-200	430 stainless steel lab-jacks	430 stainless steel	200*200
LJ550-100	304 stainless steel lab-jacks	304 stainless steel	100*100
LJ550-150	304 stainless steel lab-jacks	304 stainless steel	150*150
LJ550-200	304 stainless steel lab-jacks	304 stainless steel	200*200



BALL TYPE INTERFACE CLAMP

Item.	Name	Model	Size(mm)
CS610-12	Ball type interface clamp	12#	--
CS610-18	Ball type interface clamp	18#	--
CS610-28	Ball type interface clamp	28#	--
CS610-35	Ball type interface clamp	35#	--
CS610-50	Ball type interface clamp	50#	--



STAINLESS STEEL TAP SHAPED INTERFACE CLIP

Item.	Name	Model	Size(mm)
CC200-10	Stainless steel tap shaped interface clip	10#	--
CC200-12	Stainless steel tap shaped interface clip	12#	--
CC200-14	Stainless steel tap shaped interface clip	14#	--
CC200-19	Stainless steel tap shaped interface clip	19#	--
CC200-24	Stainless steel tap shaped interface clip	24#	--
CC200-29	Stainless steel tap shaped interface clip	29#	--
CC200-34	Stainless steel tap shaped interface clip	34#	--
CC200-40	Stainless steel tap shaped interface clip	40#	--
CC200-45	Stainless steel tap shaped interface clip	45#	--



STAINLESS STEEL SHEET INTERFACE CLAMP

Item.	Name	Model	Size(mm)
CC220-14	Stainless steel sheet interface clamp	14#	--
CC220-19	Stainless steel sheet interface clamp	19#	--
CC220-24	Stainless steel sheet interface clamp	24#	--
CC220-29	Stainless steel sheet interface clamp	29#	--



SPRING WATER STOPPER

Item.	Name	Size(mm)	Length(mm)
TC100-14	Spring water stopper	45	45
TC100-24	Spring water stopper	60	60



HUO TYPE WATER STOP CLIP (COPPER PLATING)

Item.	Name	Model	Specification
TH200-12	Huo type water stop clip (copper plating)	12#	Max.opening 12
TH200-17	Huo type water stop clip (copper plating)	17#	Max.opening 17
TH200-20	Huo type water stop clip (copper plating)	20#	Max.opening 20
TH200-25	Huo type water stop clip (copper plating)	25#	Max.opening 25
TH200-30	Huo type water stop clip (copper plating)	30#	Max.opening 30
TH200-40	Huo type water stop clip (copper plating)	40#	Max.opening 40
TH200-50	Huo type water stop clip (copper plating)	50#	Max.opening 50

SPIRAL SEALING CLIP

Item.	Name	Size(mm)	Rod(mm)
TH300-20	Spiral sealing clip (full copper)	--	--
TH400-20	The spiral sealing clip (plating)	--	--



TH300-20



TH400-20

MORSE TONGS (COPPER PLATING)

Item.	Name	Size(mm)	Rod(mm)
TC200-40	Morse tongs (copper plating)	40	--
TC200-50	Morse tongs (copper plating)	50	--
TC200-60	Morse tongs (copper plating)	60	--
TC200-70	Morse tongs (copper plating)	70	--
TC200-80	Morse tongs (copper plating)	80	--
TC201-50	Morse tongs (copper plating)	50	--
TC201-60	Morse tongs (copper plating)	60	--
TC201-70	Morse tongs (copper plating)	70	--
TC201-80	Morse tongs (copper plating)	80	--



TC200-40/TC200-50/TC200-60/TC200-70/TC200-80



TC201-50/TC201-60/TC201-70/TC201-80




GLASS TUBE CUTTER

Item.	Name	Specification	Size(mm)
CB600-30	Glass tube cutter	30	--



**BURNER**

Item.	Name	Specification	Size(mm)
BY110-10	American bunsen burner	Natural gas	--
BY110-20	American bunsen burner	Propane	--
BY110-30	American bunsen burner	Coal gas	--
BY140-10	Micro burner	Natural gas	--
BY140-20	Micro burner	Propane	--
BY140-30	Micro burner	Coal gas	--

**TRIPOD STAND**

Item.	Name	Specification	Height(mm)	Side Length(mm)
TJ100-10	Tripod stand	Paint	200	100
TJ100-125	Triangle tripod	125mm chrome finish	200	125
TJ100-150	Triangle tripod	150mm chrome finish	200	150
TJ100-200	Triangle tripod	200mm chrome finish	150	200
TJ110-125	Triangle tripod	125mm stainless steel	200	125
TJ110-150	Triangle tripod	150mm stainless steel	200	150
TJ110-200	Triangle tripod	200mm stainless steel	150	200
TJ120-125	Circular tripod	125mm chrome finish	200	125
TJ120-150	Circular tripod	150mm chrome finish	200	150
TJ120-200	Circular tripod	200mm chrome finish	150	200
TJ130-125	Circular tripod	125mm stainless steel	200	125
TJ130-150	Circular tripod	150mm stainless steel	200	150
TJ130-200	Circular tripod	200mm stainless steel	150	200

**IRON RING PAINT**

Item.	Name	Specification	Size(mm)
AA300-3*1	Iron three ring paint	--	--
AA300-60	Iron small ring paint	--	60
AA300-80	Iron middle ring paint	--	80
AA300-100	Iron big ring paint	--	100

**WITH THE TOP WIRE RING CLOSED**

Item.	Name	Size(mm)	Ring(mm)
H0330-60	With the top wire ring closed	--	60
H0330-80	With the top wire ring closed	--	80
H0330-100	With the top wire ring closed	--	100
H0330-120	With the top wire ring closed	--	120
H0330-140	With the top wire ring closed	--	140

**TOP WIRE OPENING RING**

Item.	Name	Size(mm)	Ring(mm)
H0320-50	Top wire opening ring	--	50
H0320-60	Top wire opening ring	--	60
H0320-80	Top wire opening ring	--	80
H0320-100	Top wire opening ring	--	100
H0320-120	Top wire opening ring	--	120
H0320-140	Top wire opening ring	--	140

**THE SUPPORT RING (CLOSED WITH ROD)**

Item.	Name	Ring(mm)	Length(mm)
H0310-60	The support ring (closed with rod)	60	180
H0310-80	The support ring (closed with rod)	80	195
H0310-100	The support ring (closed with rod)	100	200
H0310-120	The support ring (closed with rod)	120	215
H0310-140	The support ring (closed with rod)	140	225

**SUPPORT RING (WITH BAR OPENING)**

Item.	Name	Ring(mm)	Length(mm)
H0300-60	Support ring (with bar opening)	60	180
H0300-80	Support ring (with bar opening)	80	195
H0300-100	Support ring (with bar opening)	100	200
H0300-120	Support ring (with bar opening)	120	215
H0300-140	Support ring (with bar opening)	140	225

**STAINLESS STEEL SPOON**

Item.	Name	Specification	Head	Length
SS300-125	Stainless steel spoon	125mm	Single head	125mm
SS300-140	Stainless steel spoon	140mm	Single head	140mm
SS300-160	Stainless steel spoon	160mm	Single head	160mm
SS300-180A	Stainless steel spoon	160mm	Double head	160mm
SS300-180	Stainless steel spoon	180mm	Double head	180mm
SS300-180A	Stainless steel spoon	180mm	Double head	180mm
SS300-200	Stainless steel spoon	200mm	Single head	200mm
SS300-200A	Stainless steel spoon	200mm	Double head	200mm
SS300-220	Stainless steel spoon	220mm	Single head	220mm
SS300-220A	Stainless steel spoon	220mm	Double head	220mm
SS300-260	Stainless steel spoon	260mm	Single head	260mm

**STAINLESS STEEL DOUBLE SPATULAS, WITH TWO FLAT ENDS**

Item.	Name	Specification	Length
SP110-125	Stainless steel scraper	One flat One bent 125mm 304	125mm
SP110-150	Stainless steel scraper	One flat One bent 150mm 304	150mm
SP110-180	Stainless steel scraper	One flat One bent 180mm 304	180mm
SP110-200	Stainless steel scraper	One flat One bent 200mm 304	200mm
SP110-220	Stainless steel scraper	One flat One bent 220mm 304	220mm

**STAINLESS STEEL MICRO SPOONS, WITH SPOON AND SPATULA BLADE**

Item.	Name	Specification	Length
SP220-3*1	Stainless steel micro medicine spoon	3*1	--

STAINLESS STEEL SPOON

Item.	Name	Specification	Length
SS430-3*1	Stainless steel spoon	3*1	--
SP120-150	Stainless steel spoon	A scoop and a shovel 150mm	150mm
SP120-175	Stainless steel spoon	A scoop and a shovel 175mm	175mm
SP120-200	Stainless steel spoon	A scoop and a shovel 200mm	200mm
SP120-230	Stainless steel spoon	A scoop and a shovel 230mm	230mm
SP125-125	Stainless steel spoon	A scoop and a shovel S S 304	125mm
SP125-150	Stainless steel spoon	A scoop and a shovel S S 304	150mm
SP125-210	Stainless steel spoon	A scoop and a shovel S S 304	210mm



SS430-3*1



SP120 series



SP125 series

STAINLESS STEEL FORCEPS

Item.	Name	Specification	Length
FC200-125	Stainless steel tweezers	125mm	125mm
FC200-140	Stainless steel tweezers	140mm	140mm
FC200-160	Stainless steel tweezers	160mm	160mm
FC200-180	Stainless steel tweezers	180mm	180mm
FC200-200	Stainless steel tweezers	200mm	200mm
FC200-250	Stainless steel tweezers	250mm	250mm
FC200-300	Stainless steel tweezers	300mm	300mm



STAINLESS STEEL SCISSORS

Item.	Name	Specification	Length
SE125-140	Stainless steel scissors	140mm	140mm
SE125-160	Stainless steel scissors	160mm	160mm
SE125-180	Stainless steel scissors	180mm	180mm



STAINLESS STEEL SPATULAS

Item.	Name	Specification	Length
SP250-100	Stainless steel scraper	80mm	100mm
SP250-150	Stainless steel scraper	150mm	150mm
SP250-200	Stainless steel scraper	210mm	200mm



STERILIZING BOXES FOR PETRI DISHES

Item.	Name	Specification	Length
SC400-10	Culture dish disinfection barrel	105mm*260mm	105*260mm



SPOON

Item.	Name	Specification	Size
PP400-3*1	Plastic spoon	3*1	12.5cm/13.5cm/15cm
SP128-160	Plastic spoon	PP	160mm
SP128-180	Plastic spoon	PP	180mm
SP128-200	Plastic spoon	PP	200mm
SS400-220	Horn spoon	220mm	220mm
SS400-3*1	Horn spoon	3*1	10cm/11cm/13cm
SP100-20-1	Stainless steel burning spoon	Stainless steel	--
SP100-20-2	Copper burning spoon	Copper	--



PP400-3*1



SP128 series



SS400-220



SS400-3*1



SP100-20-1




SP100-20-2

LOOP-CARRIER / INOCULATION LOOP / VACCINATION NEEDLE / STAINLESS STEEL NEEDLE

Item.	Name	Specification	Length
KH200-10	Loop-carrier	200mm	200mm
HT100-10	Inoculation loop	50mm	50mm
HT100-20	Vaccination needle	50mm	50mm
DN-150	Stainless steel needle	stainless steel	--



KH200-10



HT100-20




DN-150

FLASK HOLDER


Item.	Name	Specification	Dia.
FB500-90	Black flask holder	Rubber	90mm
FB520-90	Red flask holder	Rubber	90mm
FB510-90	White flask holder	PP 90mm	90mm
FB510-160	White flask holder	PP 160mm	160mm
CR450-80	Round bottom flask holder	80mm	80mm
CR450-100	Round bottom flask holder	100mm	100mm
CR450-110	Round bottom flask holder	110mm	110mm
CR450-120	Round bottom flask holder	120mm	120mm
CR450-140	Round bottom flask holder	140mm	140mm
CR450-160	Round bottom flask holder	160mm	160mm





CLAY NET

Item.	Name	Specification	Size
CP850-125	Clay net	Asbestos free	125*125mm
CP850-150	Clay net	Asbestos free	150*150mm
CP850-200	Clay net	Asbestos free	200*200mm
CP850-250	Clay net	Asbestos free	250*250mm



WIRE GAUZES WITH ASBESTOS

Item.	Name	Specification	Size
AN800-125	Gauze with asbestos	Asbestos free	125*125mm
AN800-150	Gauze with asbestos	Asbestos free	150*150mm
AN800-200	Gauze with asbestos	Asbestos free	200*200mm
AN800-250	Gauze with asbestos	Asbestos free	250*250mm



DOUBLE BALL WITH NET

Item.	Name	Specification	Length
PF600-20	Double ball with net	--	--



EAR WASHING BULB

Item.	Name	Specification	Max Vol.
PF300-30	Ear Washing Bulb	--	30ml
PF300-60	Ear Washing Bulb	--	60ml
PF300-90	Ear Washing Bulb	--	90ml



HEAVY RING

Item.	Name	L.D.	Flasks
WR300-48	Heavy ring	48mm	125-700ml
WR300-56	Heavy ring	56mm	250-1000ml
WR300-61	Heavy ring	61mm	125-700ml
WR300-70	Heavy ring	70mm	700-2000ml
WR300-74	Heavy ring	74mm	1-4L



PROTECTIVE GLOVES

Item.	Name	Specification	Length
HP190-100	Protective gloves	Silica gel big size	--
HP190-90	Protective gloves	Natural rubber big size	--
HP190-80	Protective gloves	Silica gel small	--
HP190-70	Protective gloves	Natural rubber small size	--

