

CONTENTS

COMPANY PROFILE



MICROPIPETTE

P Series Whole Autoclavable Micropipette	03	Bottle Top Dispenser	07
M Series Half Autoclavable Micropipette	04	Large Capacity Pipette Pump	08
Micro Fixed Pipette	05	Pipette Pump	08
Multichannel Pipette	06	Electric Pipette	09



LAB STIRRER

Mini Magnetic Stirrer	11	Multiple Magnetic Stirrer	15
Mini Magnetic Stirrer With Display	12	Overhead Stirrer	16
Multiple Magnetic Stirrer	13	Magnetic Stirrer With Heating	18
Large Capacity Magnetic Stirrer	14	Digital Display Heating Magnetic Stirrer	19



HEATING MANTLE

Heating Mantle With Electronic Temperature Control	21	Electronic Temperature Control Heating Mantle With Magnetic stirrer	23
Digital Heating Mantle	22	Digital Temperature Control Heating Mantle With Magnetic stirrer	24



SHAKER&MIXER

Vortex Mixer	25	Multifunctional Mixer Adapted Module	29
Vortex Mixer with multifunction	27		
Digital Display Multifunctional Mixer	28	Orbital Shaker	30



LAB CENTRIFUGE

Low Speed Centrifuge	32	Digital Mini Centrifuge	34
Mini Centrifuge	33		




VACUUM PUMP

Oilless Diaphragm Vacuum Pump	36	Large Size Vacuum Filtration Apparatus	38
Oilless Diaphragm Vacuum Pump 30L	37	Vacuum Filtration Apparatus	39



DRY BATH

Mini Dry Bath	41	High Temperature Dry Bath Incubator	45
Dry Bath Incubator	43		



WATER BATH

Water Bath	48
------------	----



BALANCE SCALE

Square Pan Balance	50	0.001g Balance	52
Round Pan Balance	51	0.0001g Balance	53



INFRARED STERILIZER

Infrared Sterilizer	55
---------------------	----

P Series Whole Autoclavable Micropipette



121 °C sterilization



Wide range



Dispense accurately



TECHNICAL DATA

Model	Capacity Range	Test Capacity	Tip	Increment	Accuracy Error	Repeat Error
P10	0.5-10ul	5	10ul	0.1ul	±1.00%	±1.00%
P20	2-20ul	10	200ul	0.5ul	±1.00%	±1.00%
P50	5-50ul	25	200ul	0.5ul	±0.90%	±0.60%
P100	10-100ul	50	200ul	1ul	±1.00%	±0.40%
P200	20-200ul	100	200ul	1ul	±0.80%	±0.30%
P1000	100-1000ul	500	1000ul	5ul	±0.70%	±0.25%
P5000	1000-5000ul	2500	5000ul	50ul	±0.60%	±0.30%
P10000	2-10ml	5ml	10ml	0.1ml	±1.20%	±0.30%




M Series Half Autoclavable Micropipette

Half autoclavable



Wide range



Dispense accurately



TECHNICAL DATA

Model	Capacity Range	Test Capacity	Tip	Increment	Accuracy Error	Repeat Error
M2.5	0.1-0.25ul	1.25	10ul	0.05ul	±3.00%	±3.00%
M10	0.5-10ul	5	10ul	0.1ul	±1.50%	±0.80%
M20	2-20ul	10	200ul	0.5ul	±1.20%	±0.40%
M50	5-50ul	25	200ul	0.5ul	±0.90%	±0.30%
M100	10-100ul	50	200ul	1ul	±1.00%	±0.15%
M200	20-200ul	100	200ul	1ul	±0.70%	±0.15%
M1000	100-1000ul	500	1000ul	5ul	±0.80%	±0.20%
M5000	1000-5000ul	2500	5000ul	50ul	±0.60%	±0.15%
M10000	2-10ml	5ml	10ml	0.1ml	±1.20%	±0.20%

Micro Fixed Pipette



Smaller size and more convenient



Wide range



Dispense accurately



TECHNICAL DATA

Model	Capacity Range	Color	Test Capacity	Accuracy Error	Repeat Error	Tips
A5001	4ul	white	4ul	±1.00%	±1.00%	10ul
A5002	5ul	red	5ul	±1.20%	±1.00%	10ul
A5003	10ul	green	10ul	±1.00%	±0.80%	10ul
A5004	20ul	black	20ul	±0.60%	±0.40%	200ul
A5005	25ul	white	25ul	±0.50%	±0.40%	200ul
A5006	30ul	yellow	30ul	±0.50%	±0.40%	200ul
A5007	40ul	brown	40ul	±0.50%	±0.40%	200ul
A5008	50ul	blue	50ul	±0.50%	±0.40%	200ul
A5009	75ul	purple	75ul	±0.50%	±0.40%	200ul
A5010	100ul	orange	100ul	±0.50%	±0.40%	200ul
A5011	120ul	white	120ul	±0.50%	±0.40%	200ul



Multichannel Pipette


Digital window



8/12 channels optional



Suitable for 96-well plates



Suitable for 96-well plates



TECHNICAL DATA

Model	Type	Capacity Range	Tip	Increment	Test Capacity	Allowed maximum error
M8-10	8 channels	0.5-10ul	10ul	0.1ul	5	2.50% 0.125
M8-50	8 channels	5-50ul	200ul	0.5ul	25	1.50% 0.375
M8-300	8 channels	50-300ul	300ul	5ul	150	1.00% 1.500
M12-10	12 channels	0.5-10ul	10ul	0.1ul	5	2.50% 0.125
M12-50	12 channels	5-50ul	200ul	0.5ul	25	1.50% 0.375
M12-300	12 channels	50-300ul	300ul	5ul	150	1.00% 1.500

Bottle Top Dispenser



Digital window



Gear adjustment



Acid and alkali resistant



Safety valve protection



Various models




TECHNICAL DATA

Model	Capacity Range	Test Capacity	Tip	Increment	Accuracy Error	Repeat Error
DA-2ML	0.4-2ml	0.05	1ml	10ul	0.20%	2ul
DA-5ML	1-5ml	0.1	2.5ml	25ul	0.20%	5ul
DA-10ML	2-10ml	0.25	5ml	50ul	0.20%	10ul
DA-30ML	5-30ml	0.5	15ml	125ul	0.20%	25ul
DA-60ML	10-60ml	1	30ml	250ul	0.20%	50ul

Large Capacity Pipette Pump



Continuous inhalation



Wide range



TECHNICAL DATA

Product name

Large capacity pipette pump

Pipetting range

0.1 - 100 ml

Color

Green part & gray spherical hand pinch



Pipette Pump

Dispense easily



TECHNICAL DATA

Product name

Pipette pump

Pipetting range

2ml

10ml

25ml

Color

Blue & White

Green & White

Red & White

Electric Pipette



High quality pipette pump



Adjustable pipetting speed



Dynamic drainage
Gravity drainage




Durable power and can be
used continuously for 8 hours

TECHNICAL DATA

Product Name	Electric Pipette Pump
Range	0.1 - 100ml
Weight	190g
Operating and charging temperature	0°Cto 35°C
Pipetting speed	Pipette 50ml of liquid in 10 seconds
Battery	Lithium battery 2.4V 700mAh

LAB STIRRER



Mini Magnetic Stirrer



Strong magnetic field



Large mixing ability



Durable brushless motor



Stepless speed regulation,
wide speed range


TECHNICAL DATA

Model	MS5	Mixing direction	Right
Control method	Knob	Shell protection level	IP42
Motor type	Brushless motor	Power	5W
Working plate size	135MM	Voltage	AC100-240V/DC24V
Working plate material	PVC	Allow ambient temperature	5-4°C 80%RH
Maximum stirring volume	2L(H ₂ O)	Product Size	168*160*47MM
RPM	200-1800RPM	Package dimensions	250*200*70MM
Stir bar size	20-30MM	Gross weight / net weight	0.66/0.46KG

Mini Magnetic Stirrer With Display



Strong magnetic field



LCD display



Durable brushless motor



Large capacity

TECHNICAL DATA

Model	MS5S	Mixing direction	Right
Control method	Knob	Shell protection level	IP42
Motor type	Brushless motor	Power	5W
Working plate size	135MM	Voltage	AC100-240V/DC24V
Working plate material	PVC	Allow ambient temperature	5-4°C 80%RH
Maximum stirring volume	2L(H ₂ O)	Product Size	168*160*47MM
RPM	200-2000RPM	Package dimensions	250*200*70MM
Stir bar size	20-30MM	Gross weight / net weight	0.66/0.46KG

Digital Display Hotplate Stirrer



Strong magnetic field



Anti-corrosion ceramic plate



Accurate temperature control



Over temperature protection

TECHNICAL DATA

Model	HS5S	Working plate	Anti-corrosion ceramic plate
Control method	Brushless motor	Plate temperature	310℃
Maximum stirring volume	2L(H ₂ O)	Power	400W
RPM	50-1500RPM	Temperature sensor	Inside and outside PT1000
Stir bar size	20-40MM	Voltage	ac110 / 220V
Mixing direction	Right	Product size	155*220*700MM
Shell protection level	IP42	Net weight	0.7KG

Large Capacity Magnetic Stirrer



Timing function & lock function



Two operation modes: knob + button



Forward and reverse



Large capacity

TECHNICAL DATA

Model	MS-10L	MS-20L	MS-50L
Control method	Button + knob	Button + knob	Button + knob
Display method	LCD display	LCD display	LCD display
Power	20W	30W	60W
Maximum stirring volume	10L	20L	50L
Range of RPM	50-1500rpm	50-1500rpm	50-800rpm
Speed regulation accuracy	±10rpm/min	±10rpm/min	±10rpm/min
Timing range	23h/59min/59sec	23h/59min/59sec	23h/59min/59sec
Magnet size range	30-80mm	40-100mm	60-120mm

Multiple Magnetic Stirrer



Can be set to run synchronously or asynchronously




Timing function



Stepless speed regulation,
wide speed range



Durable brushless motor



LCD display

TECHNICAL DATA

Model	MMS4Pro	MMS6Pro	MMS8Pro	RPM	50-1500RPM
Control	Asynchronous control (synchronization can be set)			Bar size	20-50MM
Number of stirring points	4	6	8	Display type	LCD liquid crystal display
Working plate size	110MM			Timing	23h/59min/59sec
Maximum stirring volume	4*500ml	6*500ml	8*500ml	Voltage	AC100-240V
Motor type	Brushless motor			Shell protection level	IP42
Mixing direction	Right			Allow ambient temperature	5-40°C 80%RH

Overhead Stirrer



Smooth start



Adjustable speed



DC permanent magnet motor



Circuit protection

TECHNICAL DATA

Model	OS-10L	OSC-10L	OSC-20L
Maximum amount of stirring(H ₂ O)	10L	10L	20L
Motor input power	72W	72W	200W
Motor output power	60W	60W	180W
Voltage	100-240V	100-240V	100-240V
Frequency	50/60Hz	50/60Hz	50/60Hz
Power	80W	80W	210W
Range of rotate speed	100-2000rpm	100-2000rpm	100-2000rpm
Rotate speed display	/	LCD	LCD
Speed display resolution	/	±1rpm	±1rpm
Speed increment/decrement	20rpm	20rpm	20rpm
Maximum torque	20N.cm	20N.cm	40N.cm
Maximum viscosity	10000mPas	10000mPas	30000mPas
Clamping diameter range of drill chuck	0.5-10mm	0.5-10mm	0.5-10mm
Overall dimension (DxWxH)	200*315*600mm	200*315*600mm	200*315*600mm
Weight	9kg	9kg	9kg
DIN EN60529 protection mode	IP42	IP42	IP42
Allowable ambient temperature	5-40°C	5-40°C	5-40°C
Allowable ambient humidity	80%	80%	80%
Time setting	0-999min	0-999min	0-999min

Magnetic Stirrer With Heating



Strong magnetic



Aluminum plate conducts heat more uniformly


TECHNICAL DATA

Model	HS-12	HS-17	HS-19
Working plate size	12*12CM	17*17CM	19*19CM
Working plate material	Aluminum plate	Aluminum plate	Anti-corrosion ceramic plate
RPM	0-1600rpm	0-1600rpm	0-1600rpm
Function	Heating and stirring	Heating and stirring	Heating and stirring
Heating power	180W	500W	600W
Working plate temperature	Room temperature to 280°C	Room temperature to 300°C	Room temperature to 350°C
Stirring capacity	1L	3L	5L

Digital Display Heating Magnetic Stirrer



Strong magnetic



LCD display



Anti-corrosion ceramic plate



Safety circuit over temperature protection

TECHNICAL DATA

Model	Plate size	Plate material	RPM	Function	Heating power	Plate temperature	Stirring capacity
HSC-19T	19 * 19 cm	Anti-corrosion ceramic plate	200-2000 (r.p.m)	Heating and stirring	600w	Room temperature to 350 C±10%	5L

HEATING MANTLE



Digital Heating Mantle



Digital heating



Small stroke




Accurate temperature control



Heating evenly, heating up fast

TECHNICAL DATA

Model	HMC								
Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
Function	Heating								
Temperature control	Digital temperature control								
Heating power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
Max temperature	Room temperature-380 °C								



TECHNICAL DATA

Model	HMS								
Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
Function	Heating/magnetic stirrer								
Heating power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
Max temperature	Room temperature-380 ℃								



Magnetic stirring function



Heating evenly, heating up fast

TECHNICAL DATA

Model	HMSC								
Capacity	100ml	250ml	500ml	1L	2L	3L	5L	10L	20L
Function	Heating/magnetic stirrer								
Temperature control	Digital temperature control								
Heating power	100W	150W	250W	350W	450W	600W	800W	1200W	2400W
Max temperature	Room temperature-380 °C								

SHAKER&MIXER



Mini Vortex Mixer



Stepless speed regulation
Wide speed range



High quality silicone
Multicolor optional



Tap to start Compact body



Brushless motor Durable

TECHNICAL DATA

Model	VM-300S	VM-300	Shell protection level	IP42	
Control mode	Speed control knob	Fixed speed	Power	12W	
Motor	Brushless motor		Voltage	100-240V,50/60HZ	
Amplitude	4.5MM (circular oscillation)		Product Size	100*93*72MM	93*93*72MM
Maximum sample processing capacity	50ML		Allowable ambient temperature	5-40℃,80%RH	
Speed range	0-3000RPM	3000RPM	Gross / net weight	0.56/0.37KG	0.55/0.36KG

Vortex Mixer With Multifunction



Ultra-low body stable operation



Various fixtures
large handling capacity



Brushless motor with strong power



Touch/Run easy to switch

TECHNICAL DATA

Model	VM-500S	Speed range	Touch: 0-3000RPM Run: 200-3000RPM
Operating mode	Touch / Run	Shell protection level	IP21
Control mode	knob speed adjustment	Voltage	AC100-240V/DC24V
Motor	Brushless motor	Power	20W
Amplitude	4.5MM (circular oscillation)	Temperature	5-40°C, 80%RH
Max capacity	50ml*4/15ml*9/5-10ml*12/2ml*24		

Digital Display Multifunctional Mixer



LCD Display



Various fixtures, large handling capacity



Brushless motor with strong power



Touch/Run, easy to switch

TECHNICAL DATA

Model	VM-500pro	Speed range	Touch: 0-3000RPM Run: 100-3000RPM
Operating mode	Touch / Run	Timing / Time	1 - 99min / 1 - 59sec
Control mode	PWM speed adjustment	Shell protection level	IP21
Display	LCD	Voltage	AC100 - 240V / DC24V
Motor	Brushless motor	Power	20W
Amplitude	4.5MM(circular oscillation)	Temperature	5-40 °C ,80%RH
Max capacity	50ml*4/15ml*9/5-10ml*12/2ml*24		

Multifunctional Mixer Adapted Module

Model	Picture	Image	Module application scope
VM501			Centrifuge tube / Test tube
VM502			Flask / Multi tube
VM503			Φ12mm Blood collection tube * 21
VM504			15ml Centrifuge tube * 12 5ml Centrifuge tube * 12
VM505			50ml Centrifuge tube * 4 15ml Centrifuge tube * 5
VM506			0.2ml Centrifuge tube * 96 96 well plate
VM507			1.5ml Centrifuge tube * 20 2.0ml Centrifuge tube * 30
VN508			Cell culture plate

Orbital Shaker



Timing function



LCD display



Adjustable speed



Brushless motor Durable

TECHNICAL DATA

Model	OS-20	Maximum load	2KG
Speed	50-300RPM	Working temperature	4-45°C
Orbit diameter	20MM	Dimensions	284*264*90MM
Timing range	1min-99h59min	Voltage	AC100-230V

CENTRIFUGE




Low Speed Centrifuge



Easy to use



LCD display



Automatic fault alarm



Electronic door lock



TECHNICAL DATA

Model	LC-6S	LC-8S	LC-12S
Max rotation speed	500-5000r/min	500-5000r/min	500-4000r/min
Max relative centrifugal force	3074xg	3074xg	2220xg
Max capacity	6*15ml centrifuge tube	8*15ml centrifuge tube	12*5ml/10ml blood collection tube
Speed regulation accuracy	±30r/min		
Timing range	10s~99min59s		
Noise	≤60dB(A)		
Voltage	AC220V/110V-50Hz 5A		
Power	130W		
Dimensions	375*320*240 mm(L*W*H)		
Weight	5.5kg		

Mini Centrifuge



Stop running after opening lid



Multi-rotor compatible



Low noise




High frequency and wide voltage

TECHNICAL DATA

Model	MC-4S	MC-7S	MC-4/7S
Control mode	PI power direct drive		
Max rotation speed	1080xg	3286xg	1080/3286xg
Max rotation speed	4000rpm±5%	7000rpm±5%	4000/7000rpm±5%
Rotor capacity	Angled rotor 0.2/0.5/1.5/2ml*12; PCR rotor 0.2ml*12*4		
Max speed acceleration /reduction	≤12 s/ ≤15s		
Noise under max rotation speed	≤47db (A)		
Input power	20W	25W	25W
Shell protection level	AC100-240V / DC24V		
Allow surrounding temperature	5-40℃ / 80%RH		

Digital Mini Centrifuge



LCD display



Multi-rotor compatible



Rotor with big capacity



Speed adjustable and timing

TECHNICAL DATA

Model	MC-7pro	MC-12pro	MC-14.5pro
Control mode	PWM electronic speed regulation		
Display	LCD liquid crystal		
Max relative centrifugal force	3286xg	9660xg	14100xg
Max rotation speed	7000rpm±5%	12000rpm±5%	14500rpm±5%
Rotor capacity	Angled rotor 0.2/0.5/1.5/2ml*12 ; PCR rotor 0.2ml*12*4		
Timing	1-99min/1-59sec		
Max speed acceleration /reduction	≤12s/≤15s		
Noise under max rotation speed	≤47db (A)	≤56db (A)	≤56db (A)
Input power	25W	45W	60W

VACUUM PUMP



Oilless Diaphragm Vacuum Pump



Oil-free and pollution-free



Overload,overheat and short circuit protection



High precision vacuum gauge



Illuminated switch



TECHNICAL DATA

Model	VP-10L	VP-15L
Type	Negative pressure	
Suction flow rate	10L/min	15L/min
Limit Vacuum Degree	0.075Mpa	0.08Mpa
Noise	<40db	<45db
Enclosure Ratings	IP21	
Power supply	20W	30W
Diaphragm	Rubber	
Protection devices	Overload protection, over temperature protection, short circuit protection	
Gross/net weight	1.95/1.58KG	

Oilless Diaphragm Vacuum Pump - 30L



TECHNICAL DATA

Model	VP-30L
Type	Positive and negative pressure
Suction flow rate	30L/min
Limit Vacuum Degree	0.085Mpa
Greatest pressure	5.0bar
Noise	<50 (db)
Enclosure Ratings	IP21
Power supply	50(W)
Diaphragm	Rubber
Protection devices	Overload protection, over temperature protection, short circuit protection

Large Size Vacuum Filtration Apparatus



Suitable for bottles with GL-45



Large Diameter·Fast suction

TECHNICAL DATA

Receiving flask	Sand core diameter	Filter membrane	Sand core aperture
2000ml	80mm	90 - 100um	10um
5000ml	80mm	90 - 100um	10um
10000ml	80mm	90 - 100um	10um

Vacuum Filtration Apparatus



**Boro
3.3**

Boro 3.3



Tight sealing



TECHNICAL DATA

Receiving flask	Filter hole	Hole size of filter element
250ml	20mm	10um
500ml	40mm	10um
1000ml	40mm	10um
2000ml	40mm	10um
5000ml	40mm	30-50um
10000ml	40mm	30-50um


DRY BATH




Mini Dry Bath



Small volume



Multi-point operation



LCD display




Overtemperature protection

TECHNICAL DATA

Model	MBD100	MBD100-C
Temperature arrange	room temp+5°C~100°C	0°C~100°C
Temperature stability	≤±0.3°C	
Temperature Homogeneity	≤±0.3°C	
Display accuracy	0.1°C	
Heating rate	≤12min	≤12min
Cooling rate	-	≤20min
Multi-point operation	support	
Temperature point can be set	2	
Time setting range	1-999min or 1-999sec	
Power Supply	DC24V,2A	DC12V,5A
Size	162*107*110mm	162*107*126mm
Weight	0.8kg	0.9kg

Module options	Numbers and aperture of test tube
MBD1	40*0.2ml (Centrifuge tube)
MBD2	24*0.5ml (Centrifuge tube)
MBD3	15*1.5ml (Centrifuge tube)
MBD4	15*2.0ml (Centrifuge tube)
MBD5	4*5/15ml (Centrifuge tube)
MBD6	2*50ml (Centrifuge tube)
MBD7	8*12.5*12.5mm (Special module for cuvette)
OEM of modules	Can be customized based on customer requirement

Dry Bath Incubator



LCD Display



Temperature deviation calibration



Overtemperature protection



Multipoint temperature operation

TECHNICAL DATA

Model	BD100	BD100-C
Temperature arrange	Room temp+5°C~100°C	-10°C~100°C
Temperature stability	≤±0.3°C	
Module temperature Homogeneity	≤±0.3°C	
Display accuracy	0.1°C	
Heating rate	≤12min	
Cooling rate	-	≤15min (20°C to -5°C)
		≤15min (20°C to -5°C)
		Test at ambient temperature of 25°C
Multi-point operation	Support	
Temperature point can be set	5	
Time setting range	1min-99h59min	
Power Supply	AC220/110V,50/60HZ,150W	

Module options	Numbers and aperture of test tube
BD1	54*0.5ml (Centrifuge tube)
BD2	35*1.5ml (Centrifuge tube)
BD3	35*2.0ml (Centrifuge tube)
BD4	96-well plates
BD5	12*5/15ml (Centrifuge tube)
BD6	6*50ml (Centrifuge tube)
OEM of modules	Can be customized based on customer requirement

High Temperature Dry Bath Incubator



High precise temperature control



Overtemperature protection



LCD display



Temperature deviation calibration

TECHNICAL DATA

Model	BD150-1	BD150-2	BD150-4
Temperature arrange	room temp+5°C~150°C	room temp+5°C~150°C	room temp+5°C~120°C
Temperature stability (40°C)	$\leq \pm 0.3^{\circ}\text{C}$		
Temperature stability (100°C)	$\leq \pm 0.5^{\circ}\text{C}$		
Temperature Homogeneity (40°C)	$\leq \pm 0.3^{\circ}\text{C}$		
Temperature Homogeneity (100°C)	$\leq \pm 0.5^{\circ}\text{C}$		
Temperature Homogeneity	$\leq 10\text{min}$		
Display accuracy	0.1°C		
Standard equipped modules	1	2	4
Time setting range	1min-99h59min		
Power Supply	AC220/110V,200W		

Module options


Aperture

HBD1	48*0.2ml(Centrifuge tube)
HBD2	42*0.5ml(Centrifuge tube)
HBD3	24*1.5ml(Centrifuge tube)
HBD4	24*2.0ml(Centrifuge tube)
HBD5	12*5ml/15ml(Centrifuge tube)
HBD6	4*50ml(Centrifuge tube)
HBD7	96-well plates(Use two modules in combination)

WATER BATH



Water Bath



Digital display thermostatic



High precise temperature control



Stainless steel a molding tank




High corrosion resistance

TECHNICAL DATA

Model	BHS-1	BHS-2	BHS-4	BHS-6
Temperature range	Room temp-99.9℃			
Power	300W	600W	800W	1200W
Accuracy	0.1℃			
Working chamber size	16*16*13CM	30*15*14CM	30*30*14CM	45*30*13CM

BALANCE SCALE



Square Pan Balance



Unit conversion



Counting

TECHNICAL DATA

Model	JNB50001	JNB100001	JNB150001	JNB200001
Maximum capacity	5000g	10000g	15000g	20000g
Accuracy	0.1g	0.1g	0.1g	0.1g
Pan size	170*190mm			
Dimensions	260*265*75mm			

0.001g Balance



 0.001g accuracy

 Unit conversion

 Level 4 conversion

 PCS Counting

TECHNICAL DATA

Model	JNB1003G	JNB2003G	JNB3003G	JNB5003G
Maximum capacity	100g	200g	300g	500g
Accuracy	0.001g	0.001g	0.001g	0.001g
Pan size	90mm			
Dimensions	200*265*140mm			

0.0001g Balance



Internal calibration



LED display



Full electromagnetic interference protection



Multi-device printing
Standard RS232 interface

TECHNICAL DATA

Model	JAB1004 JAB1004N	JAB1204 JAB1204N	JAB2004 JAB2004N	JAB2204 JAB2204N
Maximum capacity	100g	120g	200	220g
Minimum capacity	10mg	10mg	10mg	10mg
Accuracy	0.0001g			
Pan size	80mm			

* JAB-N series is a fully automatic calibrated balance with built-in weight

Round Pan Balance



Unit conversion



Counting

TECHNICAL DATA

Model	Maximum capacity	Accuracy	Pan size	Dimensions
JNB2001	200g	0.1g	130mm	200*265*75mm
JNB5001	500g	0.1g	130mm	200*265*75mm
JNB10001	1000g	0.1g	130mm	200*265*75mm
JNB20001	2000g	0.1g	150mm	200*265*75mm
JNB30001	3000g	0.1g	150mm	200*265*75mm
JNB1002	100g	0.01g	130mm	200*265*140mm
JNB5002	500g	0.01g	130mm	200*265*140mm
JNB10002	1000g	0.01g	130mm	200*265*140mm
JNB20002	2000g	0.01g	150mm	200*265*140mm
JNB30002	3000g	0.01g	150mm	200*265*140mm
JNB1003	100g	0.001g	110mm	200*265*140mm
JNB2003	200g	0.001g	110mm	200*265*140mm

INFRARED STERILIZER



Infrared Sterilizer



Sterilization only takes 5-7 second



Can be used in an anaerobic environmen



No open-flame ,more safer



Easy to clean

TECHNICAL DATA

Model	DS-900	DS-900S
Power	150W	240W
Max. temperature inside the heating hole	825°C±50°C	
Heating temperature	Two temperatures available 400°C/800°C	
Max-diameter of sterilizer goods	Φ14mm	Φ35mm
Length of heating zone	150mm	100mm
Warm-up time	10min	

NOTE



JOANLAB®

