



BDS400 Inverted Biological Microscope

- New research and new development
- Apply to observation and cultivation for cell tissue
- Wide field eyepiece, view field up to $\Phi 23\text{mm}$, more comfortable for observation



BDS400 Inverted Biological Microscope



Long working distance condenser N.A. 0.30, Working distance: 72mm

Working distance: 195mm (without condenser), available for extra high culture dish

Large size stage, convenient for research.

Size: 240mm (X) × 210 (Y) mm

Mechanical stage available for 96 holes plate.

Moving range: 128mm (X) × 80 (Y) mm



BDS400 Inverted Fluorescence Microscope





Large diameter quintuple nosepiece can be installed more objectives, more convenient for using.

Light distribution (both): 100 : 0 (100% for eyepiece);

80 : 20 (80% for trinocular head and 20% for eyepiece)

Culture Dish Holder




Knob can be changed by left or right

Objective Parameter List (Cover-glass thickness 1.1mm)


	Objectives	N.A.	Working distance(mm)
LWD Infinity Plan Objectives	LPlanFL 4×	0.11	12.1
	LPlanFL 10×	0.25	10.3
	LPlanFL 20×	0.45	5.8
	LPlanFL 40×	0.65	5.1
LWD Infinity Plan Phase Contrast Objectives	LPlanFL PHP 4×	0.10	9.2
	LPlanFL PHP 10×	0.25	10.3
	LPlanFL PHP 20×	0.45	5.8
	LPlanFL PHP 40×	0.65	5.1




BDS400 Inverted Biological Microscope Size(mm)



BDS400 Inverted Fluorescence Microscope Size(mm)



BDS400 Inverted Biological Microscope General Layout Diagram



BDS400 Inverted Biological Microscope Outfits

Item	Specification	Outfits
Eyepiece	WF10× / 22mm	●●
LWD Infinity Plan Objectives	LPlanFL 4× / 0.11 W.D.=12.1mm	○
	LPlanFL 10× / 0.25 W.D.=10.3mm	●
	LPlanFL 20× / 0.45 W.D.=5.8mm	●
	LPlanFL 40× / 0.65 W.D.=5.1mm	●
LWD Infinity Plan Phase Contrast Objectives	LPlanFL PHP 4× / 0.10 W.D.=9.2mm	○
	LPlanFL PHP 10× / 0.25 W.D.=10.3mm	○
	LPlanFL PHP 20× / 0.45 W.D.=5.8mm	●
	LPlanFL PHP 40× / 0.65 W.D.=5.1mm	○
Annular Spot	10×/20×/ 40×	●
	4×	○
Seidentopf Trinocular Head	Inclined 45°, interpupillary distance: 48-76mm Light distribution (both): 100: 0 (100% for eyepiece) 80:20 (80% for trinocular head, and 20% for eyepiece)	●
Seidentopf Trinocular Head	Inclined 45°, interpupillary distance: 48-76mm Light distribution (both): 100: 0 and 0:100 (100% for eyepiece or 100% for trinocular head)	○
Nosepiece	Quintuple	●
Mechanical Stage	Stage size: X×Y: 210×241mm, Round slide size: Φ110mm, Attached mechanical stage (available for 96 holes plate, moving range X×Y: 128×80mm.)	●
Culture Dish Holder	65mm	○
	54mm	●
	35mm	○
	blood cell counter	○
Condenser	Long working distance, Quickly detachable, N.A.0.3, Working distance: 72mm (with condenser),195mm(without condenser).	●
Koehler Illumination	6V/30W (input voltage: 100V~240V)	●
	5WLED, input voltage: 100V~240V	○
Filter	Blue	●
	Green	●
	Amber /Grey	○
C-mount	1 X /0.5X /0.75X C-Mount (focus adjustable)	○
Epi Fluorescence Illumination	Epi fluorescence media unit, field diaphragm, center adjustable.	○
	100W mercury lamp ,5WLED(input voltage:100V~240V)	○
	B, G, V, UV fluorescence filters can be chosen	○

Note: "●"In Table Is Standard outfits, "○" Is Optional Accessories.

Packing Size: 660mm×590mm×325mm

Gross Weight: 18 kgs

Net Weight: 13.5 kgs