



超声波清洗机
Ultrasonic Bath

手机清洗注油机
Handpiece Lubrication System



WELCOME TO THE P&T TECHNOLOGY MEDICAL EQUIPMENT

 founded in 2004, is a subsidiary company of Nanyang Group, which is a high-tech modern and professional factory engaging in developing sterilizer center, having eight branch offices all over the country, to serve hospital, operating room, dentistry, ophthalmology and biomedical research institutes.

Our company is focus on Automatic Sterilizer and the relative auxiliary equipments (eg. Handpiece Lubrication and Cleaning Machine, Sealing Machine of the Sterilization Pouches, Ultrasonic Bath). All the products have passed the CE certificate via TÜV Rheinland.

Since foundation, with excellent technical talents, mature market and distribution channels and good brand reputation, we have provided our extraordinary products and services all round the world, for example Germany, France, USA, Japan, UK, Australia etc. We are looking forward to growing with our partners in together.



佰泰科技 · 非凡体验 · 成就未来


前身为宁波佰泰医疗设备有限公司，成立于2004年，并于2015年正式搬迁至浙江省衢州市常山县，是一家集科、工、贸为一体，专业从事生产、研发全自动灭菌器的高科技公司。服务于各医疗单位手术室、口腔科、眼科及生物医学科研单位。

公司主营全自动灭菌器及灭菌中心全套设备方案（手机清洗注油机、灭菌袋封口机、超声波清洗机、蒸馏水机等）。公司及产品均已通过德国TUV ISO13485和CE认证。公司研制开发的高温高压真空灭菌器采用欧洲灭菌技术（B级），通过专业工业设计与精密生产，不但使灭菌器性能及效能上一个台阶、更符合操作习惯，而且品控优异、服务有保障，广泛受到欢迎和好评。


自公司创办以来，我们本着诚信经营的原则，专业热忱的服务宗旨，勇于开拓市场。经过不懈努力拼搏，公司目前拥有优秀的技术人才、成熟的市场渠道和良好的企业信誉。我们期待进一步挖掘潜在商机，与更多合作伙伴共同成长。




CFDA Medical Device Production License
医疗器械生产许可证




CFDA Sterilizer Register Certificate
中华人民共和国灭菌器医疗器械注册证




TÜV Rheinland ISO 13485:2016
医疗器械ISO13485管理体系证书




TÜV Rheinland CE Sterilizer
Certificate
灭菌器CE证书




Certificate For Exportation of Medical
Products
中华人民共和国医疗器械出口销售证明



U.S. FDA Registration & Device Listing (Ultrasonic Cleaner)
超声波清洗机美国FDA注册列名认证



Intertek ETL Ultrasonic Bath Certificate
超声波清洗机美国ETL安全认证



Ultrasonic Bath Patent
超声波清洗机专利证书



Ultrasonic Bath (Recessed) Patent
嵌入式超声波清洗机专利证书

ENTERPRISE CULTURE

企业荣誉，再创辉煌

Talent is our most valuable resource.

We strive to attract,nurture and retain the best talent."To be Top brand" is the mission of P&T Medical since it's established.

Our reputation is built on the highest standards of professional integrity and ethic, which is our biggest asset. We never compromise on integrity.

“做中国人自己的品牌”是佰泰公司成立之初的便肩负的使命。我们的声誉建立在高尚的职业诚信和道德标准之上，这也是我们宝贵的资产之一，在此方面我们从不妥协。

压力蒸汽灭菌器 Class B Steam Sterilizer

KEY FEATURES:

- Top-Designed Reservoir
- Built-in Water Quality Sensor
- Excellent Sterilization Chamber & Rotatable Rack with 3 Trays
- Safety Automatic Electronic Lock and High Quality Door
- With Quick Release Coupling
- LCD display
- With the Function of Automatic Clean
- With Double Testing Function
- Built-in Customer Definition Program, to Expand the Scope of the Sterilized Items
- Built-in Printer / USB

产品特点：

- 国际公认蒸汽压力灭菌方式
- 方便菜单操作，用户可设置日期、时间等
- 换代式全自动电子门锁结构
- 一体成型内胆
- 高亮度液晶显示
- 内置热敏打印机，并配备有PC接口，实施远程维护
- 设置预热键，使灭菌室内温度可较快达到灭菌温度、减少灭菌程序时间
- 内置及自行可设的消毒程序，消毒物品范围广泛
- 内置净、污水存储水箱
- 可与蒸馏水系统连接自动注入蒸馏水，也可直接与污水管道连接自动排除冷凝水

功能特点：

- 欧洲B级，三次预真空
- 双重安全保护装置
- 开机自动检测
- 内置水箱
- 微电脑控制系统，液晶显示
- 7套固化灭菌程序，1套自定义程序

适用范围：

- 压力蒸汽灭菌器主要用于医疗单位手术室、口腔科、眼科及生物医学科研单位，对手术器械、齿科器械、敷料、玻璃器皿、注射器具等可耐压力蒸汽的物品进行灭菌



功能特点

上置式水箱

敞开式双水箱便于加水 and 清洗；大容量设计使一次加水可供更多次灭菌程序使用

内置水质传感器

对所加水质进行实时检测和监控，减少机器故障率，同时保护消毒器械

整体拉伸成型灭菌室和可旋转式器械架

一次性整体拉伸成型灭菌室，表面经过电抛光处理，形成一个光洁、致密的保护层。通过可旋转的托盘架，消毒腔内可使用3-5个托盘来分量消毒物品

全自动门锁和一次性压铸成型门体

换代型电子门锁结构，采用优质不锈钢铸件，对温度、压力等多重多点实时监控，设置极限保护，彻底杜绝常规手动门锁类灭菌器因操作不当导致的安全隐患

配有快接接口

位于前方的两个快速接口方便排水

高亮度液晶显示

高亮度液晶显示，美观实用。利用内置或自行设定的消毒程序，灵活灭菌，节能环保

内置自动清洗程序

自动清洗机器内管路和蒸汽发生器，避免堵塞

内置双重测试程序

内含真空度检测和泄露测试程序，让客户对机器的灭菌性能了如指掌

预置自行可设的消毒程序， 消毒物品范围广泛

内置式打印机

内置热敏打印机，直观显示灭菌过程，方便管理，配备有USB模块和PC接口实现远程维护



Trust Design Features



1

Top-Designed Reservoir

Open reservoir, convenient for water-refilling and cleaning. The reservoir is enlarged, making it available for 20 cycles use.

2

Built-in Water Quality Sensor

Inspect the quality of distilled water, reduce the failure rate and protect the sterilized instrument.

3

Excellent Sterilization Chamber & Rotatable Rack with 3 Trays

One-piece chamber is made of imported 304 stainless steel. The rotatable rack can make the sterilization chamber available for 3-5 trays.

4

Safety automatic electronic lock and high quality door

Fully automatic electronic lock. Monitoring the temperature and pressure in real time to ensure the safety

5

With Quick Release Coupling

The two quick release coupling in front of the unit make drainage easier

6

LCD display

LCD display gives the whole process of sterilization by means of words, fashion and practicalthe.

7

With the Function of Automatic Clean

Clean the inner pipe and steam generator automatically for avoiding jam after every 100 times of cycles.

8

With Double Testing Function

With Bowie & Dick testing and LEAK testing, convenient for machine performance testing.

9

Built-in Customer Definition Program, to Expand the Scope of the Sterilized Items

10

Built-in Printer

With a built-in printer, optional to record every cycle during the sterilization USB module and remote access service for option

压力蒸汽灭菌器

Class B Vacuum Steam Sterilizer

压力蒸汽灭菌器 Sterilizer BTD 23-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	2000W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 450mm
灭菌室容积	Chamber volume	23L
外形尺寸	Overall dimension (L xW xH)	600mm*475mm*430mm
净重	Net weight	60KG

压力蒸汽灭菌器 Sterilizer BTD17-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	1800W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 350mm
灭菌室容积	Chamber volume	17L
外形尺寸	Overall dimension (L xW xH)	550mm*475mm*430mm
净重	Net weight	55KG

压力蒸汽灭菌器

Class B Vacuum Steam Sterilizer

压力蒸汽灭菌器 Sterilizer BTD 12-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	1500W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	200mm * 360mm
灭菌室容积	Chamber volume	12L
外形尺寸	Overall dimension (L xW xH)	530mm*435mm*385mm
净重	Net weight	46KG

压力蒸汽灭菌器 Sterilizer BTD 8-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	1500W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	170mm * 320mm
灭菌室容积	Chamber volume	8L
外形尺寸	Overall dimension (L xW xH)	530mm*435mm*385mm
净重	Net weight	45KG

压力蒸汽灭菌器

Class B Vacuum Steam Sterilizer

压力蒸汽灭菌器 Sterilizer BTD23-A (敞开式水箱 Top Feeding)



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	2000W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 450mm
灭菌室容积	Chamber volume	23L
外形尺寸	Overall dimension (L xW xH)	640mm*475mm*430mm
净重	Net weight	58KG

压力蒸汽灭菌器 Sterilizer BTD17-A (敞开式水箱 Top Feeding)



额定工作电源	Nominal voltage	AC220 , 50HZ
额定输入功率	Nominal power	1800W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 350mm
灭菌室容积	Chamber volume	17L
外形尺寸	Overall dimension (L xW xH)	640mm*475mm*430mm
净重	Net weight	58KG

压力蒸汽灭菌器

Class B Vacuum Steam Sterilizer

手动门压力蒸汽灭菌器 Sterilizer BTS23-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	2000W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 450mm
灭菌室容积	Chamber volume	23L
外形尺寸	Overall dimension (L xW xH)	640mm*475mm*430mm
净重	Net weight	52KG

手动门压力蒸汽灭菌器 Sterilizer BTS17-A



额定工作电源	Nominal voltage	AC100/220v , 50HZ /60HZ
额定输入功率	Nominal power	1800W
额定工作温度	Sterilization temperature (user-defination)	121 °C / 134 °C (用户自定义105 °C - 134°C)
工作环境温度	Operating environment temperature	5 °C ~ 40 °C
额定工作压力	Rated working pressure	0.11 M Pa / 0.21 M Pa
灭菌室 (直径*深度)	Chamber dimension	250mm * 350mm
灭菌室容积	Chamber volume	17L
外形尺寸	Overall dimension (L xW xH)	640mm*475mm*430mm
净重	Net weight	50KG

超声波清洗机 Ultrasonic Bath

KEY FEATURES :


- User-friendly operating panel and display for time and heating
- Efficient 40 kHz ultrasonic high-performance robust transducer
- Scratch resistant, durable and easy to clean all stainless steel construction
- Degas and Sweep mode
- Overheating protection
- Cassette cleaning friendly
- Easy and fast draining of the cleaning liquid via valve and solenoid valve
- Ergonomically design recessed model

工作原理 :

- 通过换能器将功率超声频声能转换成机械振动，使清洗液中的微气泡在声波的作用下保持振动。当声压或声强受到压力到达一定程度时，气泡迅速膨胀又突然闭合的瞬间会产生冲击波，使气泡周围产生1012-1013Pa的压力及局调温，这种超声波空化产生的巨大压力能破坏不溶性污物而使他们分化于溶液中。
- “空化”现象：当超声波振动在液体中传播的音波压强达到一个大气压时，峰值可达到真空或负压，此时液体中产生一个很大的力将液体分子拉裂成空洞——空化核。空洞接近真空，在超声波压强反向达到最大时破裂产生强烈冲击将物体表面污物撞击下来。这种无数细小的空化气泡破裂产生的冲击波现象称为“空化”现象。

适用范围 :

- 本机适用于医院·医药、学校、科研、石油、化工、轻工、冶金、机械、交通、国防等行业，对仪器、仪表、电子器件、电板、半导体硅片、磁性材料、电镀件、五金件、光学镜片及零配件、音频磁头、涤纶过滤芯、喷丝板、乳胶模具、医疗器械、玻璃器皿、珠宝首饰、钟表零件、精密五金件、轴承、油嘴油泵及机械制造中零部件和几何形状复杂部件的清洗。



超声波清洗机 Ultrasonic Bath

超声波清洗机 BTX600- 3



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	200W (加热Heating) / 120W (不加热Non-Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	120 (W)
清洗槽尺寸	Cleaning tank dimension	240mm (L)*140mm (W)*100mm (H)
容积	Tank capacity	3L / 0.8gal
外形尺寸	Overall size	280mm (L)*180mm (W)*225mm (H)
净重	Net weight	4.7KG

超声波清洗机 BTX600 - 5 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	200W (加热Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	120 (W)
清洗槽尺寸	Cleaning tank dimension	240mm (L)*140mm (W)*150mm (H)
容积	Tank capacity	5L / 1.32gal
外形尺寸	Overall size	280mm (L)*180mm (W)*270mm (H)
净重	Net weight	4.5KG


超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 5** 不加热 Non-Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	120W
加热温度范围	Temperature range	/
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	120 (W)
清洗槽尺寸	Cleaning tank dimension	240mm (L)*140mm (W)*150mm (H)
容积	Tank capacity	5L / 1.32gal
外形尺寸	Overall size	280mm (L)*180mm (W)*270mm (H)
净重	Net weight	4.5KG

超声波清洗机 **BTX600 - 5** 嵌入式



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	200W (加热Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	120 (W)
清洗槽尺寸	Cleaning tank dimension	240mm (L)*140mm (W)*150mm (H)
容积	Tank capacity	5L / 1.32gal
外形尺寸	Overall size	520mm (L)*390mm (W)*360mm (H)
净重	Net weight	8.5KG

超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 7** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	350W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	220 (W)
清洗槽尺寸	Cleaning tank dimension	500mm (L)*140mm (W)*100mm (H)
容积	Tank capacity	7L / 1.85gal
外形尺寸	Overall size	530mm (L)*165mm (W)*275mm (H)
净重	Net weight	8KG

超声波清洗机 **BTX600 - 10** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	350W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	220 (W)
清洗槽尺寸	Cleaning tank dimension	300mm (L)*240mm (W)*150mm (H)
容积	Tank capacity	10L / 2.64gal
外形尺寸	Overall size	345mm (L)*285mm (W)*295mm (H)
净重	Net weight	8.5KG


超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 10** 不加热 Non-Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	220W
加热温度范围	Temperature range	/
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	220 (W)
清洗槽尺寸	Cleaning tank dimension	300mm (L)*240mm (W)*150mm (H)
容积	Tank capacity	10L / 2.64gal
外形尺寸	Overall size	345mm (L)*285mm (W)*295mm (H)
净重	Net weight	8.5KG

超声波清洗机 **BTX600 - 10** 嵌入式 Recessed



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	350W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	220 (W)
清洗槽尺寸	Cleaning tank dimension	300mm (L)*240mm (W)*150mm (H)
容积	Tank capacity	10L / 2.64gal
外形尺寸	Overall size	458mm (L)*398mm (W)*290mm (H)
净重	Net weight	12KG

超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 13** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	450W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	300 (W)
清洗槽尺寸	Cleaning tank dimension	375mm (L)*348mm (W)*150mm (H)
容积	Tank capacity	13L / 3.43gal
外形尺寸	Overall size	375mm (L)*350mm (W)*320mm (H)
净重	Net weight	11.5KG

超声波清洗机 **BTX600 - 13** 嵌入式 Recessed



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	450W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	300 (W)
清洗槽尺寸	Cleaning tank dimension	375mm (L)*348mm (W)*150mm (H)
容积	Tank capacity	13L / 3.43gal
外形尺寸	Overall size	485mm (L)*459mm (W)*330mm (H)
净重	Net weight	15.5KG

超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 19** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	450W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	300 (W)
清洗槽尺寸	Cleaning tank dimension	375mm (L)*348mm (W)*200mm (H)
容积	Tank capacity	19L / 4.94gal
外形尺寸	Overall size	375mm (L)*350mm (W)*370mm (H)
净重	Net weight	12.5KG

超声波清洗机 **BTX600 - 19** 嵌入式 Recessed



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	450W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	300 (W)
清洗槽尺寸	Cleaning tank dimension	375mm (L)*348mm (W)*200mm (H)
容积	Tank capacity	19L / 4.94gal
外形尺寸	Overall size	485mm (L)*459mm (W)*381mm (H)
净重	Net weight	16KG

超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 20** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	650W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	430 (W)
清洗槽尺寸	Cleaning tank dimension	552mm (L)*348mm (W)*150mm (H)
容积	Tank capacity	20L / 5.28gal
外形尺寸	Overall size	555mm (L)*355mm (W)*320mm (H)
净重	Net weight	16KG

超声波清洗机 **BTX600 - 20** 嵌入式 Recessed



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	650W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	430 (W)
清洗槽尺寸	Cleaning tank dimension	552mm (L)*348mm (W)*150mm (H)
容积	Tank capacity	20L / 5.28gal
外形尺寸	Overall size	665mm (L)*460mm (W)*330mm (H)
净重	Net weight	21KG

超声波清洗机 Ultrasonic Bath

超声波清洗机 **BTX600 - 40** 加热 Heating



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	1240W (加热 Heating)
加热温度范围	Temperature range	20 °C ~ 80 °C
加热定时范围	Timer setting (Mins)	1 ~ 99 (Mins)
超声波功率值	Sonic power	800 (W)
清洗槽尺寸	Cleaning tank dimension	620mm (L)*500mm (W)*150mm (H)
容积	Tank capacity	40L / 10.56gal
外形尺寸	Overall size	650mm (L)*530mm (W)*425mm (H)
净重	Net weight	25KG

手机清洗注油机

Handpiece Cleaning & Lubrication System

BTY 700



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ	
额定输入功率	Nominal power	25W	
空气压力	Working Pressure	0.30KPa(3.5-6.0kgf/cm2) (50-80psi)	
储油桶容量	Tank capacity	350ml	
外形尺寸	Overall size	285mm (L)*355mm (W)*275mm (H)	
净重	Net weight	7.5KG	
模式	手机类型	清洁注油时间	空气清洁时间
短模式-S	手动涡轮手机	35 秒	60 秒
	使用NAKANSHI NM系列机头的手机		
	使用NAKANSHI SGM系列机头的手机		
短模式-L	弯角手机	45秒	80秒
短模式-E	其他手机（短模式及长模式不适配的手机）	50秒	90秒

功能特点：

- 可随时简单、准确地对直弯或高速手机进行清洗和注油
- 安全性高，操作简单、易于掌握、任何人都可进行操作
- 能同时给三支口腔手机清洗及注油
- BTY700的旋转齿轮能从不同角度给手机有效地清洗和注油
- 完成清洗和注油步骤后，按下。空气键。能去除手机内多余残
- 可选择短模式、长模式、超长模式对不同手机进行注油
- 油雾过滤棉能防止油雾的扩散减少到最小程度

蒸馏水机 Water Distiller

HR-01



额定工作电源	Nominal voltage	AC100/220 , 50HZ / 60HZ
额定输入功率	Nominal power	750 W
温度保护器	Temperature Protector	160 °C
容量	Tank Capacity	4L
蒸馏量	Distillation Capacity	1L/H
净重	Net Weight	4KG

功能特点：

- 蒸馏水机被广泛应用于医院、牙科诊所、健康中心合格机构研究机构
- 蒸馏水机也可应用于灭菌器
- 请注意纯净水并不等同于蒸馏水