

## FUNCTION

- 5-inch color LCD display, color graphic menu, man-machine dialogue, convenient and intuitive
- 12 leads synchronous acquisition, 12 channel ECG waveforms in same screen display and simultaneous printing
- Manual and automatic printing of standard 12 channels, 3channels + rhythm leads, 6 channels, 6 channel + rhythm leads, and 60 second arrhythmia analysis
- Large capacity ECG data storage, support SD memory card capacity expansion
- The electrode falls off alarm and indicates the electrode falling off part
- V1, V3, V5 chest physical examination mode, easy to operate, very suitable for large-scale physical examination
- Eliminates baseline drift and other disturbances with a unique high-precision digital filter, does not cause distortion to the ECG waveform, greatly enhances the anti-baseline drift capability, and facilitates waveform interpretation
- Store 250 ECG data files



Button Operation



Quick and easy operation, save time

Data Storage(Optional)



Large data storage, support for capacity expansion

**KX-1203**  
ECG Monitor

