

FUNCTION

- 7-inch color LCD display, color graphic menu, man-machine dialogue, convenient and intuitive
- 12 leads synchronous acquisition, 12 channel ECG waveforms in same screen display and simultaneous printing
- Touch screen operation, quick and easy operation, save time
- Manual and automatic printing of standard 12 channels, 3channels + rhythm leads, 6 channels, 6 channel + rhythm leads, and 60 second arrhythmia analysis
- Large capacity ECG data storage, support SD memory card capacity expansion
- The electrode falls off alarm and indicates the electrode falling off part
- V1, V3, V5 chest physical examination mode, easy to operate, very suitable for large-scale physical examination
- Eliminates baseline drift and other disturbances with a unique high-precision digital filter, does not cause distortion to the ECG waveform, greatly enhances the anti-baseline drift capability, and facilitates waveform interpretation
- Store 250 ECG data files



Button Operation



Quick and easy operation, save time

Data Storage(Optional)



Large data storage, support for capacity expansion

Touch Screen



Quick and easy operation, save time

KX-1206
ECG Monitor

