

## FUNCTION



- ◆ 10.1-inch color LCD display, color graphic menu, man-machine dialogue, convenient and intuitive
- ◆ 12 leads synchronous acquisition, 12 channel ECG waveforms in same screen display and simultaneous printing
- ◆ Touch screen operation, quick and easy operation, save time
- ◆ Manual and automatic printing of standard 12 channels, 3 channels + rhythm leads, 6 channels, 6 channel + rhythm leads, and 60 second arrhythmia analysis
- ◆ Large capacity ECG data storage, support SD memory card capacity expansion
- ◆ The electrode falls off alarm and indicates the electrode falling off part
- ◆ V1, V3, V5 chest physical examination mode, easy to operate, very suitable for large-scale physical examination
- ◆ Eliminates baseline drift and other disturbances with a unique high-precision digital filter, does not cause distortion to the ECG waveform, greatly enhances the anti-baseline drift capability, and facilitates waveform interpretation
- ◆ Store 250 ECG data files

## PROFESSION

Button Operation




Quick and easy operation, save time

Data Storage(Optional)



Large data storage, support for capacity expansion

Interface figure



Color graphic menu, man-machine dialogue, convenient and intuitive

Touch Screen




Quick and easy operation, save time

Accessories



Lead Wire



Clip



Sucking ball

**KX-1212**  
ECG Monitor

