

KX8000 Multi-parameter monitor series

Quality assurance/safe and effective/customer service

KX8000 series multi-parameter monitor

Classic fusion, intelligent design

Large font display interface
Breathing oxygen graph display interface
Standard display interface
Total lead display interface
Unique remote monitor
Trend display interface
7-leads ECG display in same screen
Voice alarm
Analog output supported (optional)
Lithium batteries supported
Working duration for 150 minutes



12.1 inch TFT high resolution color display

KX-8000C multi-parameter monitor



Powerful software function

Centralized alarm allocation function
Large font display function
Optional waveform parameters display area
Analog signal output function



Bracket solution

Several flexible installation solutions, including the wall bracket, cart and locking bracket.



Abundant monitoring parameters

Basic configuration: ECG, NIBP, RESP, SpO2, TEMP and PR
Optional configuration: CO2, IBP, CO, AG, 12-
leads, touch screen



Personalized appearance and structure design

High resolution TFT, can display 12 waveforms simultaneously
The touch screen control makes the operation more quick and simple, and time and labor saving
The internal frame structure and optimized line design provides powerful waterproof and anti-drop ability and facilitates the maintenance.

Large amount of information

The host itself provides the holographic recording function.

The holographic information includes:

- All collected waveform data
- All parameter data
- All alarm information

It provides the multi-level multiple sound and light alarm function.

/// The printed picture is only for reference for there may be slight color deviation between it and the product. The color of the physical products shall prevail. ///