


# Fighting against COVID-19, we are together.

With the massive global outbreak of COVID-19, Aeonmed customized this respiratory work station to save more lives according to the clinical requirements.



## Shangrila 510S

Integrated Respiratory Work Station

With advanced technology and versatile ventilation functions, 510S respiratory work station meets the breathing requirements of adult and pediatric.



## Shangrila 510S

Integrated Respiratory Work Station



### Complete Ventilation Modes

- The device hosts six ventilation modes.
  - Besides the conventional Volume-Controlled ventilation mode, it also has Pressure-Controlled ventilation mode,
  - SIMV and CPAP modes.
- The machine can provide a precise tidal volume.
- It also provides inspiration hold and manual breath, one-button operation, and is thus incredibly convenient.

### Ideal O<sub>2</sub> Delivery

- With its advanced technology, the range of FiO<sub>2</sub> can be set from 40% to 100%.
- Furthermore, this range can be adjusted continuously.

### Clinical Enhancements

- Manual breath
- Inspiratory hold
- Alarm silence
- Pre-use test

### Working Scenarios

- Ward, Emergency Department,  
Sub-acute Care Unit