

HL-A310A

Manual-hydraulic Operation Table



Ot Table Traction Equipment Series



Led Operation Shadowless Lamp





Medical Pendant



乐孚科技
Hi-Life Technology

OPERATION
ROOM
SOLUTION



HL-B320A

Electro-hydraulic Operating Tables



MULTI-FUNCTIONAL OPERATION TABLE

Suitable for hand surgery, orthopedics, cardiothoracic surgery, general surgery, urology, obstetrics and gynecology, ophthalmology, plastic surgery, rectal surgery, otolaryngology and endoscopic surgery, etc.

TECHNICAL REQUIREMENTS

ITEM	TECHNICAL REQUIREMENTS
PLATFORM DIMENSION	Length≥2070cm Width≥52cm
HI-LO RANGE	54cm~104cm
BACK PLATE	Up≥80° Overall Lateral Titling≥30° Down≥40° Overall Trendelenburg & Reverse–Trendelenburg≥20°
HEADREST	Up≥15° Down≥90°
FOOT PLATE	Up≥15° Down≥90°
LOADING CAPACITY	≥250kg

Action:



HL-B320B

Electric Operating Table



T-shaped base

TECHNICAL REQUIREMENTS

ITEM	TECHNICAL REQUIREMENTS
LENGTH	≥2070mm
LOWEST HEIGHT	≤700mm
OVER WIDTH	≥550mm
HI-LO RANGE	≥300mm
TRENDELENBURG & REVERSE–TRENDELENBURG	≥25°
LATERAL TILTING	≥15°
KIDNEY BRIDGE ELEVATION	≥110mm
BACK PLATE	Up≥75° Down≥20°
HEADREST	Up≥45° Down≥90°
FOOT PLATE	Down≥90° Outward90° Up≥15°
HORIZONTAL SLIDING	≥300mm



Standard Type

To meet the needs of surgeries, including chest, abdomen, ENT, obstetrics and gynecology, urology and orthopedics.