

内窥镜系列手术器械

ENDOSCOPIA MEDICAL INSTRUMENT



MEDICAL INSTRUMENT

- 胸腹腔镜手术器械
- 妇科手术器械
- 耳鼻喉科手术器械
- 泌尿外科手术器械
- 骨科手术器械
- 附属设备

产品目录

公司简介 Company profile

钱塘江尽到桐庐，水碧山青画不如。

桐庐恒丰医疗器械有限公司位于秀丽的富春江畔、中国最美县桐庐。

公司主要生产腹腔外科、胸腔外科、妇产科、耳鼻喉科、泌尿外科、骨科等领域产品。在临床应用中受到广大专家、教授好评。

. locates in Tonglu, which is one of the most beautiful counties in China.

Our products cover many areas, including Abdominal Surgery, Thoracic Surgery, Obstetrics & Gynecology, Otolaryngological Surgery, Urological Surgery and Orthopedic Surgery. All of them are highly praised by professionals.

公司宗旨

珍惜至高的生命

奉献完美的产品

精益求精，精细工作，精做客户满意产品及优质的售后服务

The purpose of our company is to appreciate precious lives and offer perfect products.

In the meantime, higher quality products and better after-sale service are always we chasing after.

公司信念

客户的需要是我们努力的方向

客户的满意是我们追求的目标

愿与各位医疗界同仁共同努力，致力于医疗器械 “桐庐特色，中国制造”

We believe that the need of customer is our developing direction and the satisfaction of customer is our objective.

We long for cooperate with other companies to put effort on "characteristics of Tonglu, made in China"

目录 Contents

胸腹腔镜手术器械 LAPAROSCOPE

腹腔镜、穿刺器 Laparoscope, Trocar	01
钛夹、钛夹钳、施夹钳 Titanium Clip Applier	02-03
可转弯手术器械 Flexible Instrument	03
单孔（经脐部）腹腔镜手术器械 Single Hole Instruments (Umbilical cord)	04
手术钳、剪刀 Forceps and Scissors	04-09
切开刀 Bile Duct Knife	09
穿刺针（胆管针）、转换器、气腹针 Needle, Reducer, Veress Needle	10
疝气针、穿刺缝合针 Bile Duct Needle, Protective Obturator Blunt	10
双极电凝钳、双极分离钳 Double-electrode Forceps	11
单极电凝线、接口 Single Eletrode Cable	11
腹腔镜甲状腺手术器械 Thyroid Operation Instruments	12
金手指拉钩、撑开器（二叶）、（三叶）、荷包抓钳 Forceps, Retractor, Purse Stitching Clamp	12
弹簧抓钳、胆囊造影钳、可拆换吸引器 Spring Grasping Forceps, Bile Duct Forceps, Suction&Irrigation	13
电凝、带点凝吸引器 Monopolar HF Electrode, Monopolar HF Electrode with S&I	14
持针钳 Needle-holder Forceps	14
胸腔镜成套手术器械（双关节）Thoracoscopic Instruments	15
胆道取石抓钳、悬吊式腹腔镜手术器械 Biliary Tract-Lithotomy Forceps, Suspensory Laparoscopic Instruments	16
三叶扇形钳、五叶扇形钳 Fan-shaped Forceps	17
可弯取物钳 Esophagus Retractor	17
腹腔镜弯、直肠手术器械 Colon and Rectum Instruments	17
一次性使用医疗器械 Disposable Instruments	18

妇科手术器械 GYNAECOLOGY

HF6518F型 电动子宫切除器及配套手术器械 Electric Uterus-cutter and Uccessories	19
HF6518-I 电动子宫切除器及手术器械 Electric Uterus-cutter and Accessories	19
举宫器 Multi-functional Uterine Manipulator Set	20-21
HF6518F.8 阴式子宫切除手术器械（非脱垂子宫经阴道手术用）Hysteroscopy Set Apparatus	22
专用阴道缝合拉钩 Special Vaginal Retractor(Suture Needle)	22
宫腔内窥镜、宫腔活检钳、宫腔镜 Hysterosocopy Set, Uterine Biopsy Forceps, Hysteroscopy Set	23

泌尿外科手术器械 UROLOGICAL INSTRUMENTS

HF6518M-SSJ 膀胱碎石镜 Bladder Stone Endoscope	24
HF6518M-PGJ 尿道膀胱镜（插口式）Urethro-cystoscopy Set (Click)	24
HF6518M-DGJ型 前列腺汽化电切镜 Vaporization Resectoscope	25
输尿管肾镜 Ureterorenoscopy Set	25

耳鼻喉科手术器械 E.N.T INSTRUMENTS

耳内窥镜、中耳钳 Otoendoscope, Middle Ear Forceps 70mm	26
耳科手术器械 Otoscopy Instruments	27-28
耳用手术电钻 Otoscope Cutting Aspirator	28
鼻窦手术钳 Nasal Operating Forceps	29-30
鼻窦刨削器（鼻电钻）Uptodate Oesophagoscope	30
鼻科手术器械 Nasal Instruments	31
支撑喉镜 Self-retaining Laryngoscopy Set	32
显微喉钳、剪 Laryngeal Forceps, Laryngeal Scissors	33
食道镜手术器械 Esophagoscopy Instruments	34
带光纤小儿支气管镜 Pediatric Bronchoscopy Set	34

肛肠科手术器械 INTESTINE INSTRUMENTS

乙状结肠镜 Sigmoidoscope	35
---------------------	----


骨科手术器械 ORTHOPEDICS

关节镜手术器械 Arthroscopy Set	35
椎间孔镜手术系统 Transforaminal Endoscopic Surgical System	36

附属设备 MEDICAL EQUIPMENTS

内窥镜接口 Endoscopy Interface	37
医用器械消毒盒 Medical Instruments Sterilizing Box	37
电子、电气设备 Endoscopic Equipments	38-39
模拟训练箱 Simulation Training Box	40
耳鼻喉诊疗仪 Therapeutice Equipment	40
治疗椅、医生椅 Medichair, Doctor-chair	40

腹腔内窥镜 Laparoscope



穿刺器 Trocar

		穿刺器 (十字密封帽) Cross type HF6518.006.1 $\Phi 5.5 \times 105\text{mm}$ HF6518.007.1 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (十字密封帽带保护) Cross type with Safety HF6518.006.2 $\Phi 5.5 \times 105\text{mm}$ HF6518.007.2 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (磁性片阀) Magnetism(flat) HF6518.008.1 $\Phi 5.5 \times 105\text{mm}$ HF6518.009.1 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (磁性片阀带螺纹) Magnetism(flat) with screw HF6518.008.3 $\Phi 5.5 \times 105\text{mm}$ HF6518.009.3 $\Phi 10.5 \times 105\text{mm}$
		转换器 (wolf式) Reducer HF6518.089.4 $\Phi 10\text{mm}-\Phi 5\text{mm}$
		密封帽 Sealing cap HF6518.100.1 $\Phi 3\text{mm}$ 、 $\Phi 5.5\text{mm}$ 、 $\Phi 10.5\text{mm}$
		密封帽 (十字) Sealing cap (Cross type) HF6518.100 $\Phi 3\text{mm}$ 、 $\Phi 5.5\text{mm}$ 、 $\Phi 10.5\text{mm}$
		密封圈 Sealing ring HF6518.100.2 $\Phi 5.5\text{mm}$ 、 $\Phi 10.5\text{mm}$
		穿刺器 (磁性球阀) Magnetism (ball) 6518.012.1 $\Phi 5.5 \times 105\text{mm}$ 6518.013.1 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (磁性球阀带保护) Magnetism (ball) with Safety 6518.012.2 $\Phi 5.5 \times 105\text{mm}$ 6518.013.2 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (按钮式) Button type HF6518.015.1 $\Phi 5.5 \times 105\text{mm}$ HF6518.015.2 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (按钮式) Button type HF6518.015.1 $\Phi 5.5 \times 105\text{mm}$ HF6518.015.2 $\Phi 10.5 \times 105\text{mm}$
		穿刺器 (翻盖式) Multifunctional valve HF6518.010.1 $\Phi 5.5 \times 105\text{mm}$ HF6518.011.1 $\Phi 10.5 \times 105\text{mm}$

钛夹、钛夹钳 Titanium Clip Applier



钛夹
HF6518Y.575 大号4#

钛夹
HF6518Y.576 中扁3#

钛夹
HF6518Y.577 中号2#

钛夹
6518Y.578 小号1#



钛夹钳 Titanium Clip Applier
HF6518.001A $\Phi 10 \times 330\text{mm}$ (单动 大号)



钛夹钳 Titanium Clip Applier
HF6518.001B $\Phi 10 \times 330\text{mm}$ (单动 中号)



钛夹钳 Titanium Clip Applier
HF6518.001C $\Phi 10 \times 330\text{mm}$ (单动 中扁)



钛夹钳 Titanium Clip Applier
HF6518.002 $\Phi 10 \times 330\text{mm}$ (双动 中扁)



钛夹钳 Titanium Clip Applier
HF6518.003 $\Phi 10 \times 330\text{mm}$ (直角 中扁)



钛夹钳 Titanium Clip Applier
HF6518.004 $\Phi 5 \times 330\text{mm}$ (单动 小号)

施夹钳 (开放式) Clip Applier (Open type)



HF6518.005A $\Phi 10 \times 330\text{mm}$
(配塑料夹 金色) 加大号



HF6518.005B $\Phi 10 \times 330\text{mm}$
(配塑料夹 紫色) 大号



HF6518.005C $\Phi 5 \times 330\text{mm}$
(配塑料夹 绿色) 中大号



施夹钳 Clip Applier



施夹钳 Tyco Clip Applier (配泰科可吸收结扎夹)
HF6518.005 $\Phi 10 \times 330\text{mm}$

配泰科可吸收结扎夹
Tyco Ligation Clip

施夹钳、结扎夹 Clip Applier , Polymer Ligation



结扎夹（塑料夹）加大号
HF6518Y.570 金色



结扎夹（塑料夹）大号
HF6518Y.571 紫色



结扎夹（塑料夹）中大号
HF6518Y.572 绿色



施夹钳 Clip Applier
HF6518.005.2 $\Phi 10 \times 330\text{mm}$
加大号（金黄色）

施夹钳 Clip Applier
HF6518.005.3 $\Phi 10 \times 330\text{mm}$
大号（紫色）

施夹钳 Clip Applier
HF6518.005.4 $\Phi 5 \times 330\text{mm}$
中大号（绿色）

施夹钳 Clip Applier



施夹钳 Clip Applier
HF6518.005.5 $\Phi 10 \times 330\text{mm}$
（配蛇牌）



施夹钳 Clip Applier
HF6518.005.1 $\Phi 10 \times 330\text{mm}$
（配强生吸收夹）



可转弯手术器械 Flexible Instrument

分离钳（可转弯）
HF6518.802.1 $\Phi 5 \times 410\text{mm}$




直剪刀（可转弯）
HF6518.801 $\Phi 5 \times 410\text{mm}$

分离钳（可转弯）
HF6518.802 $\Phi 5 \times 410\text{mm}$

无创抓钳（可转弯）
HF6518.803 $\Phi 5 \times 410\text{mm}$

单孔（经济部）腹腔镜手术器械 Single Hole Instruments (Umbilical cord)



手术钳、剪刀 Forceps and Scissors

手术钳、剪拆卸图

拆卸步骤：按图示序号①逆方向旋转转轮；②张开手柄；③向外拉钳管和钳芯；④旋转钳芯；⑤向外拉钳芯。
安装步骤：⑤④③②①



手术钳、剪刀 Forceps and Scissors

	直剪 (单动) Straight-Scissors HF6518.021 $\Phi 5 \times 330\text{mm}$ HF6518.020 $\Phi 10 \times 330\text{mm}$		弯分离钳 (30°) Marryland forceps(30°curved) HF6518.028.5 $\Phi 5 \times 330\text{mm}$
	直剪 (双动) Straight-Scissors HF6518.021.1 $\Phi 5 \times 330\text{mm}$ HF6518.020.1 $\Phi 10 \times 330\text{mm}$		弯分离钳 (45°) Marryland forceps(45°curved) HF6518.028.6 $\Phi 5 \times 330\text{mm}$
	弯剪 (双动) Curved-Scissors HF6518.022 $\Phi 5 \times 330\text{mm}$		直角分离钳 Dissecting forceps(90°curved) HF6518.029 $\Phi 5 \times 330\text{mm}$
	弯剪 (微型) Curved-Scissors HF6518.213 $\Phi 5 \times 330\text{mm}$		肠钳 (中空) Grasping forceps (Large-jaw/Fenestrated) HF6518.037 $\Phi 5 \times 330\text{mm}$
	直剪 (上翘) Straight-Scissors HF6518.023 $\Phi 5 \times 330\text{mm}$		肠钳 (中空) Grasping forceps (Large-jaw/Fenestrated) HF6518.037.1 $\Phi 5 \times 330\text{mm}$
	钩剪 Hook-Scissors HF6518.024 $\Phi 5 \times 330\text{mm}$		肠钳 (无齿) Grasping forceps (Large-jaw/toothless) HF6518.038 $\Phi 5 \times 330\text{mm}$
	直分离钳 Dissecting forceps(straight) HF6518.027 $\Phi 5 \times 330\text{mm}$		肠钳 (有齿) Grasping forceps(Large-jaw) HF6518.038.1 $\Phi 5 \times 330\text{mm}$
	弯分离钳 Marryland forceps (curved) HF6518.028 $\Phi 5 \times 330\text{mm}$		无损伤抓钳 (中空) Atraumatic forceps(Fenestrated) HF6518.039 $\Phi 5 \times 330\text{mm}$
	弯分离钳 (大弯) Marryland forceps(curved) HF6518.028.1 $\Phi 5 \times 330\text{mm}$		鸭嘴抓钳 (加长) Grasping forceps(Duck-mouth) HF6518.039.1 $\Phi 5 \times 330\text{mm}$
	弯分离钳 Marryland forceps(curved) HF6518.028.2 $\Phi 5 \times 330\text{mm}$		无创抓钳 Atraumatic forceps HF6518.040 $\Phi 5 \times 330\text{mm}$
	弯分离钳 (微型双动) Marryland forceps(curved) HF6518.028.3 $\Phi 5 \times 330\text{mm}$		无创伤抓钳 Atraumatic forceps HF6518.040.1 $\Phi 5 \times 330\text{mm}$
	弯分离钳 (70°) Marryland forceps(70°curved) HF6518.028.4 $\Phi 5 \times 330\text{mm}$		无损伤抓钳 Atraumatic forceps HF6518.041 $\Phi 5 \times 330\text{mm}$

手术钳、剪刀 Forceps and Scissors

	无损伤抓钳 Atraumatic forceps HF6518.041.1 $\Phi 5 \times 330\text{mm}$		粗齿抓钳 (三孔) Triplicate forceps HF6518.054 $\Phi 5 \times 330\text{mm}$
	无损伤抓钳 (粗齿) Atraumatic forceps(Rough tooth) HF6518.042 $\Phi 5 \times 330\text{mm}$		活检钳 Small Biopsy forceps HF6518.056 $\Phi 5 \times 330\text{mm}$ 长头 Long head HF6518.056.1 $\Phi 5 \times 330\text{mm}$ 短头 Short head
	重型抓钳 Maxi forceps HF6518.043 $\Phi 5 \times 330\text{mm}$		带针活检钳 Biopsy forceps HF6518.058 $\Phi 5 \times 330\text{mm}$ 长头 Long head HF6518.058.1 $\Phi 5 \times 330\text{mm}$ 短头 Short head
	倒齿抓钳 Grasping forceps(Inverted tooth) HF6518.044 $\Phi 5 \times 330\text{mm}$		胆囊抓钳 Babcock grasping forceps HF6518.059 $\Phi 5 \times 330\text{mm}$
	倒齿抓钳 (加长) Grasping forceps (Inverted tooth, lengthen) HF6518.045 $\Phi 5 \times 330\text{mm}$		输卵管抓钳 Oviduct forceps HF6518.061 $\Phi 5 \times 330\text{mm}$
	鼠牙抓钳 (2+3齿 双动细齿) Rat toothed forceps HF6518.047 $\Phi 5 \times 330\text{mm}$		输尿管抓钳 (有齿) Ureteral grasping forceps HF6518.061.1 $\Phi 5 \times 330\text{mm}$
	鼠牙抓钳 (1+2齿 双动无齿) Rat toothed forceps HF6518.047.1 $\Phi 5 \times 330\text{mm}$		输尿管抓钳 (无齿) Ureteral grasping forceps HF6518.061.2 $\Phi 5 \times 330\text{mm}$
	鼠牙抓钳 (2+3齿 双动粗齿) Rat toothed forceps HF6518.048 $\Phi 5 \times 330\text{mm}$		系膜抓钳 Mesentery forceps HF6518.062 $\Phi 5 \times 330\text{mm}$
	鼠牙抓钳 (2+3齿 双动粗齿加长) Rat toothed forceps HF6518.048.2 $\Phi 5 \times 330\text{mm}$		系膜抓钳 (中空) Mesentery forceps HF6518.062.1 $\Phi 5 \times 330\text{mm}$
	鼠牙抓钳 (2+3齿 单动细齿) Rat toothed forceps HF6518.049 $\Phi 5 \times 330\text{mm}$		横齿抓钳 Grasping forceps HF518.063 $\Phi 5 \times 330\text{mm}$
	钝头抓钳 Blunt grasping forceps HF6518.050 $\Phi 5 \times 330\text{mm}$		竖齿抓钳 Grasping forceps HF6518.064 $\Phi 5 \times 330\text{mm}$
	粗齿抓钳 (中空) Grasping forceps (Rough tooth/Fenestrated) HF6518.053 $\Phi 5 \times 330\text{mm}$		O型抓钳 (中空) O-Grasping forceps(Fenestrated) HF6518.066 $\Phi 5 \times 330\text{mm}$

手术钳、剪刀 Forceps and Scissors


	阑尾抓钳 Appendiceal grasping forceps HF6518.067 $\Phi 5 \times 330\text{mm}$		细齿抓钳 (单动/2孔) Croce-olmi grasping forceps HF6518.099 $\Phi 5 \times 330\text{mm}$
	阑尾抓钳 Appendiceal grasping forceps HF6518.067.1 $\Phi 5 \times 330\text{mm}$		抓钳 (小抓钳) Grasping forceps HF6518.135 $\Phi 5 \times 330\text{mm}$
	胃抓钳 Stomach grasping forceps HF6518.071 $\Phi 5 \times 330\text{mm}$		抓钳 (有齿碗型) Grasping forceps HF6518.136 $\Phi 5 \times 330\text{mm}$
	肝抓钳 Liver grasping forceps HF6518.073 $\Phi 5 \times 330\text{mm}$		抓钳 (带钉钝头) Grasping forceps HF6518.137 $\Phi 5 \times 330\text{mm}$
	吻合牵引钳 Tally pull forceps HF6518.074 $\Phi 5 \times 330\text{mm}$		抓钳 (无损伤 长头) Grasping forceps HF6518.138 $\Phi 5 \times 330\text{mm}$
	B型打结钳 (S型) B-knot forceps HF6518.078 $\Phi 5 \times 330\text{mm}$		抓钳 (肺叶式) Grasping forceps HF6518.139 $\Phi 5 \times 330\text{mm}$



A柄 A type



B柄 B type



C柄 C type



D柄 D type



E柄 E type



F柄 F type

手术钳、剪刀 Forceps and Scissors

	大肠钳 Grasping forceps(Large-jaw) HF6518X.008 Φ5×330mm		大型直角分离钳 Dissecting forceps(90°) HF6518.026.1 Φ10×330mm
	大型抓取钳 (单动) Large grasping forceps HF6518.018.4 Φ10×330mm		大型直角分离钳 Dissecting forceps(90°) HF6518.026.2 Φ10×330mm
	单齿大抓钳 (1+1齿) Large grasping forceps (1+1 tooth) HF6518.018.3 Φ10×330mm		肠钳 Grasping forceps HF6518.038.2 Φ10×330mm
	取石钳 Biopsy forceps HF6518.019 Φ10×330mm		倒齿大抓钳 Grasping forceps (Inverted tooth) HF6518.045.1 Φ10×330mm
	直剪 (单动) Straight-Scissors HF6518.020 Φ10×330mm		胆囊大抓钳 Babcock grasping forceps HF6518.059.1 Φ10×330mm
	弯剪 (双动) Curved-Scissors HF6518.022.2 Φ10×330mm		胆囊大抓钳 Babcock grasping forceps HF6518.059.2 Φ10×330mm
	大型弯分离钳 Dissecting forceps(curved) HF6518.025 Φ10×330mm		胃大抓钳 Stomach grasping forceps HF6518.071.1 Φ10×330mm
	大型直角分离钳 Dissecting forceps(90°) HF6518.026 Φ10×330mm		胃大抓钳 Stomach grasping forceps, HF6518.071.2 Φ10×330mm

手术钳、剪刀 Forceps and Scissors

	胃圆头抓钳 Stomach grasping forceps HF6518.072 Φ10×330mm		板克抓钳Grams clamp plate HF6518.126.1 Φ10×330mm
	吻合牵引钳Tally pull forceps HF6518.074.1 Φ10×330mm		板克抓钳Grams clamp plate HF6518.126.2 Φ10×330mm
	宽嘴无齿抓钳 Grasping forceps (wide mouth toothless) HF6518.503 Φ10×330mm		缝合钳Suture Forceps HF6518.130 Φ10×330mm
	血管阻断钳 Vascular occlusion forceps HF6518.079 Φ10×330mm		卵圆钳 Oval forceps HF6518.132 Φ10×330mm

取石钳（开放式）Biopsy forceps
HF6518.019.3 Φ10×250mm



切开刀 Bile Duct Knife



胆总管切开刀 Bile duct knife
6518.080 Φ5×330mm



胆总管切开刀 Bile duct knife
6518.081 Φ10×330mm

穿刺针（胆管针） Needle



转换器 Reducer

HF6518.089 插管式 $\Phi 10-\Phi 5\text{mm}$



HF6518.089A 100mm



HF6518.089B 50mm



HF6518.089C 22mm



转换器（卡式）Reducer
HF6518.089.1 $\Phi 10-\Phi 5\text{mm}$



转换器（翻盖式）Reducer
HF6518.089.3 $\Phi 10-\Phi 5\text{mm}$

气腹针 Veress Needle

气腹针 Veress Needle

HF6518.090 $\Phi 2 \times 120\text{mm}$
HF6518.090.1 $\Phi 2.5 \times 120\text{mm}$
HF6518.090.2 $\Phi 3.0 \times 120\text{mm}$




气腹针接头 Adaptor of Veress Needle
6518.097.1（中）



气腹针接头 Adaptor of Veress Needle
6518.097.2（小）

疝气针、穿刺缝合针 Bile Duct Needle, Protective Obturator Blunt




穿刺缝合针 Suture needle for puncture
6518.120.1 疝气针 $\Phi 2.2\text{mm}$
6518.120 疝气针 $\Phi 1.5\text{mm}$



穿刺缝合针 Suture needle for puncture
6518.123 疝气针 $\Phi 1.5\text{mm}$




穿刺缝合针 Suture needle for puncture
6518.023.1 疝气针 $\Phi 1.5\text{mm}$



穿刺缝合钳 Suture forceps for puncture
6518.119 疝气钳 $\Phi 2.5\text{mm}$
6518.119.1 疝气钳 $\Phi 2.5\text{mm}$

双极电凝钳 Double-electrode Forceps



双极分离钳 Double-electrode Forceps




单极电凝线、接口 Single Eletrode Cable



腹腔镜甲状腺手术器械 Thyroid Operation Instruments



金手指拉钩 Forceps



撑开器 (二叶)、(三叶) Retractor



荷包抓钳 Purse Stitching Clamp



弹簧抓钳 Spring Grasping Forceps



弹簧抓钳 (短头)
Spring grasping forceps(Short head)
HF6518.068 $\Phi 5 \times 330\text{mm}$

弹簧抓钳 (长嘴)
Spring grasping forceps(Long head)
HF6518.069 $\Phi 5 \times 330\text{mm}$

弹簧抓钳 (长嘴中空)
Spring grasping forceps(Long head, hollow)
HF6518.069.1 $\Phi 5 \times 330\text{mm}$

胆道造影钳 Bile Duct Forceps



胆道造影钳 Bile duct forceps
HF6518.055 $\Phi 5 \times 330\text{mm}$

可拆换吸引器 Suction&Irrigation



HF6518.034 三通式 $\Phi 5/\Phi 10 \times 330\text{mm}$

HF6518.034.2 史托斯式 $\Phi 5 \times 330\text{mm}$

HF6518.034.1A 妇科专用吸引管 $\Phi 5 \times 330\text{mm}$

HF6518.034.1B 妇科专用吸引管 $\Phi 10 \times 330\text{mm}$

HF6518.035.1 盒式按钮 $\Phi 5 \times 330\text{mm}$

HF6518.035 双按钮 $\Phi 5/\Phi 10 \times 330$

HF6518.035.2 单按钮 $\Phi 5/\Phi 10 \times 330\text{mm}$

电凝 Monopolar HF Electrode



电钩 Hook-Electrode
HF6518.082 $\Phi 5 \times 330\text{mm}$



电棒 Ballpoint-Electrode
HF6518.083 $\Phi 5 \times 330\text{mm}$



电铲 Spatula-Electrode
HF6518.084 $\Phi 5 \times 330\text{mm}$



电钩(陶瓷头) Hook-Electrode
HF6518.082.1 $\Phi 5 \times 330\text{mm}$



电刀 Knife-Electrode
HF6518.086 $\Phi 5 \times 330\text{mm}$



电针 Needle-Electrode
HF6518.085 $\Phi 5 \times 330\text{mm}$

带电凝吸引器 Monopolar HF Electrode with S&I



电凝(钩)吸引器 Hook type
HF6518.088 $\Phi 5 \times 330\text{mm}$



滑动式 Driven



电凝(棒状)吸引器 Rodlike
HF6518.035.3 $\Phi 5/\Phi 10 \times 330\text{mm}$




可转弯吸引器 Turning type $\Phi 5 \times 410\text{mm}$

持针钳 Needle-holder Forceps

持针钳 V型(重型) V-Needle-holder forceps
HF6518.033A $\Phi 5 \times 330\text{mm}$



持针钳 O型 O-Needle-holder forceps
HF6518.032 $\Phi 5 \times 330\text{mm}$



持针钳(直头)
HF6518.031



持针钳(弯头)
HF6518.031.1



持针钳(自动复位)
HF6518.031.2



持针钳(自动复位夹线)
HF6518.031.3



HF6518.031持针钳(V型)
Needle-holder forceps
 $\Phi 5 \times 330\text{mm}$



HF6518.032持针钳(O型)
Needle-holder forceps
 $\Phi 5 \times 330\text{mm}$




HF6518.033C持针钳(枪型)
Needle-holder forceps
 $\Phi 5 \times 330\text{mm}$


胸腔镜成套手术器械 (双关节) Thoracoscopic Instruments



胆道取石抓钳 Biliary Tract-Lithotomy Forceps



悬吊式腹腔镜手术器械 Suspensory Laparoscopic Instruments




三叶扇形钳、五叶扇形钳 Fan-shaped Forceps



可弯取物钳 Esophagus Retractor



腹腔镜弯、直肠手术器械 Colon and Rectum Instruments



一次性使用医疗器械 Disposable Instruments



一次性套管穿刺器
HF6518Y.552 Φ10mm



HF6518Y.568 手术缝合针线



一次性套管穿刺器
HF6518Y.553 Φ12mm



一次性微创筋膜缝合器 HF6518Y.564



一次性套管穿刺器
HF6518Y.551 Φ5mm



一次性气腹针
HF6518Y.556 Φ2.1mm



一次性使用手术电极
HF6518Y.566 DJ-3
HF6518Y.567 DJ-4



一次性标本袋
HF6518Y.563



一次性取物袋
HF6518Y.565

HF6518F型 电动子宫切除器及配套手术器械 Electric Uterus-cutter and Uccessories



HF6518F.023.1 电动马达线



电机 Import drive Motor
HF6518F.023 电动马达



HF6518F.022 控制器 Control unit



穿刺套管 Trocar sleeve
HF6518F.008 $\Phi 15 \times 95\text{mm}$
HF6518F.009 $\Phi 18 \times 95\text{mm}$




转换器 Reducer
HF6518F.013 $\Phi 15-5\text{mm}$
HF6518F.014 $\Phi 18-5\text{mm}$



扩张器 Dilator
HF6518F.010 $\Phi 10/\Phi 15\text{mm}$
HF6518F.011 $\Phi 10/\Phi 18\text{mm}$



宫颈钳 Cervical forceps
HF6518F.021.1 250mm



宫颈钳 Cervical forceps
HF6518F.021 250mm



腹式切除刀管 Cutting tube (TAH)
HF6518F.001 $\Phi 10 \times 320\text{mm}$
HF6518F.002 $\Phi 15 \times 320\text{mm}$
HF6518F.003 $\Phi 18 \times 320\text{mm}$




腹式切除刀管 Cutting tube (TAH)
HF6518F.003.1 $\Phi 15 \times 320\text{mm}$



阴式切除刀管 Cutting tube(LAVH)
201.004 $\Phi 15 \times 320\text{mm}$
201.005 $\Phi 18 \times 320\text{mm}$



HF6518F.007 拔棒 Moving bar $\Phi 5\text{mm} \times 420\text{mm}$



HF6518F.006 量棒 Uterus measuring stick $\Phi 5\text{mm} \times 500\text{mm}$



HF6518F.012 引导棒 Guide probe $\Phi 10\text{mm} \times 330\text{mm}$



肌瘤钻 Hysteromyoma drill
HF6518F.019 $\Phi 5\text{mm} \times 420\text{mm}$
HF6518F.020 $\Phi 10\text{mm} \times 420\text{mm}$



HF6518F.015 推结器 Knot-pushing stick $\Phi 5\text{mm} \times 330\text{mm}$



子宫小抓钳 Grasping forceps
HF6518F.016 $\Phi 10\text{mm} \times 420\text{mm}$
HF6518F.017 $\Phi 5\text{mm} \times 420\text{mm}$



HF6518F.018 举宫器 Uterine manipulator $\Phi 5\text{mm} \times 400\text{mm}$

HF6518-I 电动子宫切除器及手术器械 Electric Uterus-cutter and Accessories



宫颈钳 Cervical forceps
HF6518F.021.1 250mm



宫颈钳 Cervical forceps
HF6518F.021 250mm



子宫小抓钳 Grasping forceps
HF6518F.016 $\Phi 10\text{mm} \times 420\text{mm}$
HF6518F.017 $\Phi 5\text{mm} \times 420\text{mm}$



转换器 Reducer
HF6518F.013 $\Phi 15\text{-}5\text{mm}$
HF6518F.014 $\Phi 18\text{-}5\text{mm}$



扩张器 Dilator
HF6518F.010 $\Phi 10/\Phi 15\text{mm}$
HF6518F.011 $\Phi 10/\Phi 18\text{mm}$



腹式切除刀管 Cutting tube (TAH)
HF6518F.101 $\Phi 18 \times 330\text{mm}$
HF6518F.102 $\Phi 15 \times 330\text{mm}$
HF6518F.103 $\Phi 10 \times 330\text{mm}$



穿刺套管 Trocar sleeve
HF6518F.109 $\Phi 18 \times 102\text{mm}$
HF6518F.110 $\Phi 15 \times 102\text{mm}$




刀管内芯 Inner core of cutting tube
HF6518F.101.1 $\Phi 18\text{mm}$
HF6518F.102.1 $\Phi 15\text{mm}$



转换杆 Reducer
HF6518F.104 $\Phi 10\text{-}5\text{mm}$
HF6518F.105 $\Phi 15\text{-}10\text{mm}$
HF6518F.106 $\Phi 18\text{-}10\text{mm}$



HF6518F.023.1 电动马达线



HF6518F.006 量棒 Uterus measuring stick $\Phi 5\text{mm} \times 500\text{mm}$



HF6518F.007 拔棒 Moving bar $\Phi 5\text{mm} \times 420\text{mm}$



HF6518F.012 引导棒 Guide probe $\Phi 10\text{mm} \times 330\text{mm}$



HF6518F.015 推结器 Knot-pushing stick $\Phi 5\text{mm} \times 330\text{mm}$



HF6518F.018 举宫器 Uterine manipulator $\Phi 5\text{mm} \times 400\text{mm}$



肌瘤钻 Hysteromyoma drill
HF6518F.019 $\Phi 5\text{mm} \times 420\text{mm}$
HF6518F.020 $\Phi 10\text{mm} \times 420\text{mm}$




HF6518F.023 电动马达
电机 Import drive Motor




HF6518F.022 控制器 Control unit


简易举宫器 Multi-functional Uterine Manipulator




多功能举宫器 Multi-functional Uterine Manipulator Set



多功能举宫器 Multi-functional Uterine Manipulator Set




多功能子宫操纵器 Multi-functional Uterus Apparatus




HF6518F.8 阴式子宫切除手术器械 (非脱垂子宫经阴道手术用) Hysterectomy Set Apparatus




专用阴道缝合拉钩 Special Vaginal Retractor(Suture Needle)




宫腔内窥镜 Hysteroscopy Set



宫腔活检钳 Uterine Biopsy Forceps




宫腔镜 Hysteroscopy Set




HF6518M-SSJ 膀胱碎石镜 Bladder Stone Endoscope




HF6518M-PGJ 尿道膀胱镜 (插口式) Urethro-cystoscopy Set (Click)




HF6518M-DGJ型 前列腺汽化电切镜 Vaporization Resectoscope



输尿管肾镜 Ureterorenoscopy Set



耳内窥镜 Otoendoscope



耳内窥镜 Otoendoscope
HF6518E.01 0° Φ3×90mm



耳内窥镜 Otoendoscope
HF6518E.02 30° Φ3×90mm



耳内窥镜 Otoendoscope
HF6518E.03 70° Φ3×90mm



中耳钳 Middle Ear Forceps 70mm



HF6518E.001 中耳息肉钳, 咬切头
(头长/5mm)
Punch head
(head length 5mm)



HF6518E.002 中耳息肉钳
咬切头杯形头 右弯 45°
Cup head, right-curved



HF6518E.003 中耳息肉钳, 杯形头
(头长, 5mm) 45°
Cup head
(head length 5mm)



中耳息肉钳, 杯形头 Cup head
HF6518E.004.1 大 (头长/5mm)
Large (head length 5mm)
HF6518E.004.2 小 (头长/4mm)
Small (head length 4mm)



HF6518E.005 中耳息肉钳, 环形头
(头长/5mm)
Ring head
(head length 5mm)



中耳息肉钳, 麦粒头 Triticeous head
HF6518E.006.1 大号 (头长/5mm)
Large (head length 5mm)
HF6518E.006.2 小号 (头长/3.5mm)
Small (head length 3.5mm)
HF6518E.006.3 上翘 upward 30°
HF6518E.006.4 头长 head length 11mm



HF6518E.007 耳息肉剪
(直/头长/4mm)
Polyp scissors, straight
(head length 4mm)



耳息肉剪 Polyp scissors
HF6518E.008.1 (左弯/头长/4mm)
left curved (head length 4mm)
HF6518E.008.2 (右弯/头长/4mm)
right curved (head length 4mm)



HF6518E.009 镫骨钳
(头长/5mm)
Stapes install forceps
(head length 5mm)



镫骨足弓剪 Malleus arch scissors
HF6518E.010.1 左上翘 left upward
HF6518E.010.2 右上翘 right upward



HF6518E.011 槌骨剪(80mm)
Malleus rongeur scissors



HF6518E.012 槌骨剪(110mm)
Malleus rongeur scissors



HF6518E.013 耳钳, 橄榄头(头长/3.5mm)
Olive head
(head length 5mm)



HF6518E.014 耳钳, 直碗型(头长/3.5mm)
Straight bowl
(head length 3.5mm)



HF6518E.015 耳钳(头长/3.5mm)
Triticeous head
(head length 3.5mm)




HF6518E.016 耳钳, 直剪型(头长/3.5mm)
Straight scissors
(head length 3.5mm)



HF6518E.036
中耳显微息肉钳

耳科手术器械 Otoscopy Instruments



耳科手术器械 Otoscopy Instruments



乳突牵开器(固定式)
Plester self retaining retractors
HF6518E.033 130mm



乳突牵开器
Plester self retaining retractors
HF6518E.035



乳突牵开器
Backman laminectomy retractors 160mm
HF6518E.034 活动式

耳用手术电钻 Otoscope Cutting Aspirator



HF6518E.102 控制器
Control unit



HF6518E.103 脚踏开关
Foot pedal

HF6518E.101 手机马达



HF6518E.104 直接头 Straight joint



HF6518E.105 弯接头 Curved joint



耳钻头 球形 Ear drill



HF6518E.106 Φ2×70mm



HF6518E.107 Φ3×70mm



HF6518E.108 Φ4×70mm




HF6518E.109 Φ5×70mm




HF6518E.110 Φ6×70mm

鼻窦手术钳 Nasal Operating Forceps







鼻窦手术钳 Nasal Operating Forceps



鼻窦刨削器 (鼻电钻) Uptodate Oesophagoscope



鼻科手术器械 Nasal Instruments



支撑喉镜 Self-retaining Laryngoscopy Set



喉镜
Laryngoscope tube
HF6518H.001.1 大(L) 175mm
HF6518H.001.2 小(S) 175mm



喉镜(配0°镜)
Laryngoscope tube
for 0° endoscope
HF6518H.002.1 大(L) 175mm
HF6518H.002.2 小(S) 175mm



喉镜, 直达
Direct laryngoscope tube
HF6518H.004.1 大(L) 16×175mm
HF6518H.004.2 中(M) 14×175mm
HF6518H.004.3 小(S) 12×175mm



HF6518H.51 70° Φ8×170mm



HF6518H.005 支撑架 Hest support system



HF6518H.15 手柄 Handle



吸引管 Suction
HF6518H.16 Φ4×25mm
HF6518H.17 Φ3×25mm
HF6518H.18 Φ2.5×25mm



HF6518H.12.1 喉
刀 (镰刀状)



HF6518H.12.2
喉刀 (棱形三角头)



HF6518H.12.3
喉刀 (圆形弯30°)



HF6518H.12.4
喉刀 (半圆垂直切割)



HF6518H.12.5
喉刀 (圆形直向切割)



扁桃体刀 一柄三头
Tonsil guillotine, three heads
HF6518H.401.1
HF6518H.401.2




喉剪 Laryngeal scissors
HF6518H.44 360°旋转



粘膜喉钳 Laryngeal forceps
HF6518H.45 360°旋转

显微喉钳、剪 Laryngeal Forceps, Laryngeal Scissors



	喉剪 Laryngeal scissors	HF6518H.01	直型
	喉剪 Laryngeal scissors	HF6518H.02	上弯45°
	喉剪 Laryngeal scissors	HF6518H.02.1	上弯30°
	喉钳 Laryngeal forceps	HF6518H.04	粘膜型、倒三角
	喉钳 Laryngeal forceps	HF6518H.04.1	左弯、倒三角
	喉钳 Laryngeal forceps	HF6518H.04.2	右弯、倒三角
	喉钳 Laryngeal forceps	HF6518H.05	直、碗形口
	喉钳 Laryngeal forceps	HF6518H.06	上翘、碗形口
	喉钳 Laryngeal forceps	HF6518H.07	左开、碗形口
	喉钳 Laryngeal forceps	HF6518H.08	右开、碗形口
	喉钳 Laryngeal forceps	HF6518H.09	三角 直型
	喉钳 Laryngeal forceps	HF6518H.10	三角 上翘
	喉剪 Laryngeal scissors	HF6518H.11	左弯
	喉剪 Laryngeal scissors	HF6518H.12	右弯
	喉钳 Laryngeal forceps	HF6518H.19	左开、卵圆口
	喉钳 Laryngeal forceps	HF6518H.20	右开、卵圆口

食道镜手术器械 Esophagoscopy Instruments



HF6518H.213 食道镜手柄
Esophagoscope-handle



食道镜 Esophagoscope
HF6518H.201 $\Phi 9 \times 13 \times 250\text{mm}$
HF6518H.202 $\Phi 11 \times 15 \times 300\text{mm}$
HF6518H.203 $\Phi 10 \times 16 \times 300\text{mm}$
HF6518H.204 $\Phi 13 \times 17 \times 380\text{mm}$
HF6518H.205 $\Phi 12 \times 18 \times 400\text{mm}$



吸引管 Suction tube
6518H.210 $\Phi 2.5 \times 450\text{mm}$
6518H.211 $\Phi 3 \times 450\text{mm}$
6518H.212 $\Phi 4 \times 450\text{mm}$



食道异物钳 Foreign body forceps
HF6518H.206 鸭口形 $\Phi 4 \times 500\text{mm}$
HF6518H.207 鸭口形 $\Phi 4 \times 350\text{mm}$



食道异物钳 Foreign body forceps
HF6518H.208 倒齿形 $\Phi 4 \times 500\text{mm}$
HF6518H.209 倒齿形 $\Phi 4 \times 350\text{mm}$



带光纤小儿支气管镜 Pediatric Bronchoscopy Set

小儿支气管镜 Pediatric bronchoscopy set
HF6518Z.101 $\Phi 7 \times 280\text{mm}$
HF6518Z.102 $\Phi 6 \times 280\text{mm}$
HF6518Z.103 $\Phi 5 \times 260\text{mm}$
HF6518Z.104 $\Phi 4 \times 250\text{mm}$
HF6518Z.105 $\Phi 3.5 \times 250\text{mm}$



吸引管 Suction tube
HF6518Z.109 $\Phi 2.5 \times 320\text{mm}$
HF6518Z.110 $\Phi 3 \times 320\text{mm}$
HF6518Z.111 $\Phi 4 \times 320\text{mm}$




气管异物钳 Trachea foreign body forceps
花生头形 Peanut shape head
HF6518Z.106 $\Phi 2 \times 360\text{mm}$


气管异物钳 Trachea foreign body forceps
麦粒形 Triticeous head
HF6518Z.107 $\Phi 2 \times 360\text{mm}$

气管异物钳 Trachea foreign body forceps
鳄鱼嘴头形 Crocodile-mouth shape head
HF6518Z.108 $\Phi 2 \times 360\text{mm}$


乙状结肠镜 Sigmoidoscope



关节镜手术器械 Arthroscopy Set



椎间孔镜手术系统 Transforaminal Endoscopic Surgical System



内窥镜接口 Endoscopy Interface 导光束 Light Guide Cable

导光束 Light guide cable

HF6518.231 普通/螺口 $\Phi 4 \times 1800\text{mm}$
HF6518.232 普通/插口 $\Phi 4 \times 1800\text{mm}$
HF6518.233 普通/插口 $\Phi 4 \times 2000\text{mm}$
HF6518.234 普通/插口 $\Phi 4 \times 2500\text{mm}$
HF6518.235 普通/插口 $\Phi 4 \times 3000\text{mm}$
HF6518.236 普通/插口 $\Phi 4 \times 3500\text{mm}$
HF6518.237 普通/插口 $\Phi 5 \times 2500\text{mm}$
HF6518.238 普通/插口 $\Phi 5 \times 3000\text{mm}$
HF6518.239 耐高温/插口 $\Phi 4 \times 2500\text{mm}$
HF6518.240 耐高温/插口 $\Phi 4 \times 3000\text{mm}$

HF6518.241 耐高温/插口 $\Phi 4 \times 3500\text{mm}$
HF6518.242 耐高温/插口 $\Phi 5 \times 2500\text{mm}$
HF6518.243 耐高温/插口 $\Phi 5 \times 3000\text{mm}$
HF6518.244 配STORZ冷光源
HF6518.245 配OLYMPUS冷光源
HF6518.246 配WOLF内窥镜
HF6518.247 配STRYKER内窥镜
HF6518.909 普通/螺口 $\Phi 4 \times 1800\text{mm}$
HF6518.909.1 普通/插口 $\Phi 5 \times 1800\text{mm}$
HF6518.910 普通/螺口 $\Phi 4 \times 2500\text{mm}$
HF6518.910.1 普通/插口 $\Phi 5 \times 2500\text{mm}$
HF6518.911 普通/插口 $\Phi 5 \times 3000\text{mm}$



医用器械消毒盒 Medical Instruments Sterilizing Box



HF6518.970 膀胱镜消毒盒
For Cystoscopy set
530×230×90mm



HF6518.973 电切镜消毒盒
For Resectoscope set
420×275×90mm



HF6518.977 鼻窦手术器械消毒盒
For Sinuscopy instruments
500×245×70mm



HF6518.971 鼻窦手术器械消毒盒
For Sinuscopy set
360×255×90mm




HF6518.974 电子宫切割器消毒盒
For electric hystera-cutter
560×275×140mm



HF6518.978 内窥镜、器械通用
Endoscope Instruments in common use
500×245×70mm



HF6518.972 腹腔镜消毒盒
Endoscope laparoscope
450×90×45mm



HF6518.975 椎间盘(孔)镜消毒盒
For discectomy endoscope
350×120×70mm



HF6518.980 内窥镜消毒盒
For endoscope
290×80×50mm



HF6518.984 腹腔镜手术器械消毒盒
For laparoscope
460×80×50mm
HF6518.981 输尿管肾镜消毒盒
For ureterorenoscopy
720×12×70mm



HF6518.976 腹腔镜手术器械消毒盒
For laparoscope
540×230×70mm



HF6518.982 耳鼻喉内窥镜消毒盒
For E.N.T endoscope
300×240×52mm

全高清腹腔镜摄像系统



数字高清影像工作站

显示器 22" LCD
DVD刻录机
内窥镜接口
计算机主机 DELL i5-2400
打印机 HP6000彩喷
脚踏开关
高档铝合金多功能台车
高清腹腔镜专用软件



HF6518.914 全高清一体式内窥镜摄像系统
Integrated Endoscopic Imaging system



HF6518.900 医用全高清摄像机
Medical Camera



HF6518.903 高频电刀
High frequency electric-cutter



HF6518.902 JRQ-1 CO₂ 气腹机

电子、电气设备 Endoscopic Equipments



HF6518.904 ZLG-1型氙灯冷光源
Cold light source for hernia lamp



HF6518.905 L100型氙灯冷光源
Cold light source for hernia lamp



HF6518.905.1 便携式冷光源
Portable cold light source



医用加压器 (冲吸泵) Endoscopy pump
HF6518.906.1 JRG-1 JC200



HF6518.906 L100 医用加压器 (冲洗泵)
Medical pressure pump



医用灌注泵 Irrigation pump
HF6518.906.2 JRG-1 Jg200



HF6518.907 医用手术器械清洗器
Surgical instrument cleaning machine



ZZY-III 医用加压器 (灌注泵)
HF6518.906.3

模拟训练箱 Simulation Training Box

腹腔镜手术训练模拟配置清单：
模拟箱一套 监视器一套 持针钳一把 抓钳一把 弯分离钳一把 剪刀一把



HF6518.908 腹腔镜手术训练模拟箱

耳鼻喉诊疗仪 Therapeutice Equipment

标准配置

主操作台	1套
压缩机马达（正压泵）	1台
真空吸引马达（负压泵）	1台
喷雾枪	3把
预热除雾装置	1套
弱电控制强电的控制器	1套
主付吸引瓶	4只
机械臂	1套
器械盘	1套
废物灌	2只
药棉灌	5只
万向照明灯	1套
吸引头	1套

可选配置

- 高亮度冷光源
- 医生专用观片灯
- 诊疗椅
- 医生椅
- 头灯
- 内窥镜及手术器械
- 内窥镜图像系统
- 手术显微镜
- 显微镜图像系统

治疗椅、医生椅 Medichair, Doctor-chair

升降 Disappearing bumper 530mm-730mm
俯仰 Pitching 90°-180°
旋转 Whirl 360°
尺寸 Size 840mm×640mm×530mm或730mm
重量 Weight 125kg



医生椅 Doctor-chair



治疗椅 Medichair

